



## ORDER FORM



To see QUICKIE Power Options  
or Accessories, scan here or visit:  
[https://www.sunrisemedical.com/  
power-wheelchair-options-accessories](https://www.sunrisemedical.com/power-wheelchair-options-accessories)




**Q300M** with SEDEO® PRO™  
and JAY® Syncra™ Seating

**EIPW40**

| Customer  |  | Ship To |     |
|-----------|--|---------|-----|
| Account # |  | Name    |     |
| Date      |  | Address |     |
| PO#       |  | City    |     |
| Buyer     |  | State   | Zip |
| Mark For  |  | Phone   |     |

## Client Measurements

|                            |  |    |
|----------------------------|--|----|
| a) Top of head to seat     |  | in |
| b) Top of shoulder to seat |  | in |
| c) Underarm to seat        |  | in |
| d) Elbow to seat           |  | in |
| e) Top of leg to seat      |  | in |
| f) Shoulder width          |  | in |
| g) Chest width             |  | in |
| i) Hip width               |  | in |
| j) Actual sitting depth    |  | in |
| m) Foot length             |  | in |
| n) Knee to foot            |  | in |
| o) Seat width              |  | in |



## Notes

## Chair Specifications

|                                          |                                                                                            |
|------------------------------------------|--------------------------------------------------------------------------------------------|
| HCPC                                     | K0848 / K0856 / K0861                                                                      |
| Weights starting at:<br>(with Batteries) | NPC: 287 lbs<br>SPC: 312 lbs<br>MPC: 322 lbs                                               |
| Power Seating                            | 50° of CG Tilt and Power Recline with<br>4" mechanical sheer reduction / 12" seat elevator |
| Seat Type                                | SEDEO® PRO™ / JAY® Syncra™                                                                 |
| Seat Size Range                          | Width 16" to 22"<br>Depth 14" to 22"                                                       |
| Seat-to-Floor Height Range               | 16.5" to 19"                                                                               |

## Mid-Wheel Drive Base Specifications

|                        |                                                                    |
|------------------------|--------------------------------------------------------------------|
| Drive Wheel Position   | Mid-wheel drive                                                    |
| User Weight Capacity   | 300 lbs.                                                           |
| Base Size              | Width 20.5", Length 33"                                            |
| Minimum Turning Radius | 17" (excludes footplates)                                          |
| Drive Control          | R-net PM 90A / VR-2 90A / LiNX 120A                                |
| Driving Range          | Up to 17 miles                                                     |
| Drive Motors           | 4-pole 6 mph / 10 kph                                              |
| Drive Wheels           | 12" w/flat-free inserts (black)<br>14" w/flat-free inserts (black) |
| Casters                | 7" front / 7" rear (black)                                         |
| Charger                | 8A Dual Mode (Compact)                                             |
| Battery Size           | 50Ah Group 22                                                      |
| Color                  | 9 shroud colors and 7 color accents                                |
| Ground Clearance       | 3"                                                                 |

| BASE                                                                                 |          |       |                                                                                                                                                                                                                                                                                                                                                                 |                                                                    |            |
|--------------------------------------------------------------------------------------|----------|-------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------|------------|
| OPTION                                                                               |          | HCPC  | DESCRIPTION                                                                                                                                                                                                                                                                                                                                                     | RETAIL                                                             |            |
| Group 3 Standard Duty Mid Wheel Drive Base <i>End Item Number EIPW40 (PRO™ Seat)</i> |          |       |                                                                                                                                                                                                                                                                                                                                                                 | \$                                                                 |            |
| <input type="checkbox"/>                                                             | 400PC291 | K0848 | Quickie Q300 M NPC No Power Option for SEDEO® PRO™ Seat<br>Included in Base Pricing : Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".     | \$8,541.00                                                         |            |
| <input type="checkbox"/>                                                             | 400PC292 | K0856 | Quickie Q300 M SPC Single Power Option for SEDEO® PRO™ Seat<br>Included in Base Pricing : Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5". | \$8,541.00                                                         |            |
| <input type="checkbox"/>                                                             | 400PC293 | K0861 | Quickie Q300 M MPC Tilt/Recline for SEDEO® PRO™ Seat<br>Included in Base Pricing: Q300 M Base, Non-Expandable Electronics, No Power Seat Functions, Six-Wheel Suspension, 12" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".         | \$8,541.00                                                         |            |
| User Weight                                                                          |          |       |                                                                                                                                                                                                                                                                                                                                                                 |                                                                    | \$         |
| <input type="checkbox"/>                                                             | 400WG23  | —     | Up to 300 lbs.                                                                                                                                                                                                                                                                                                                                                  | NC                                                                 |            |
| Motor Package                                                                        |          |       |                                                                                                                                                                                                                                                                                                                                                                 |                                                                    | \$         |
| <input type="checkbox"/>                                                             | 400PF99  | —     | 4-pole 10 kph / 6 mph                                                                                                                                                                                                                                                                                                                                           | NC                                                                 |            |
| Tracking                                                                             |          |       |                                                                                                                                                                                                                                                                                                                                                                 |                                                                    | \$         |
| <input type="checkbox"/>                                                             | 400MI01  | K0108 | R-net                                                                                                                                                                                                                                                                                                                                                           | SureTrac® Module & Motor Encoders <i>(Improves drive tracking)</i> | \$1,628.00 |
| <input type="checkbox"/>                                                             | 400MI04  | K0108 | LiNX                                                                                                                                                                                                                                                                                                                                                            | LiNX Traction Control System for LiNX <i>(GYRO Technology)</i>     | \$595.00   |
| Drive Wheels                                                                         |          |       |                                                                                                                                                                                                                                                                                                                                                                 |                                                                    | \$         |
| <input type="checkbox"/>                                                             | 400DW13  | —     | Drive Wheel, 12", Solid, Knobby, Black                                                                                                                                                                                                                                                                                                                          | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400DW58  | —     | Drive Wheel, 14", Solid, Knobby, Black                                                                                                                                                                                                                                                                                                                          | NC                                                                 |            |
| Caster Wheels                                                                        |          |       |                                                                                                                                                                                                                                                                                                                                                                 |                                                                    | \$         |
| <input type="checkbox"/>                                                             | 400CW59  | —     | Caster Wheels 7", Solid, Black                                                                                                                                                                                                                                                                                                                                  | NC                                                                 |            |
| Shroud Color                                                                         |          | NEW   | Base Shroud and Back Shroud are painted to match.                                                                                                                                                                                                                                                                                                               |                                                                    | \$         |
| <input type="checkbox"/>                                                             | 400C164  | —     | Black                                                                                                                                                                                                                                                                                                                                                           | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C165  | —     | Candy Apple Red                                                                                                                                                                                                                                                                                                                                                 | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C166  | —     | Candy Blue                                                                                                                                                                                                                                                                                                                                                      | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C167  | —     | Candy Purple                                                                                                                                                                                                                                                                                                                                                    | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C176  | —     | Green Apple                                                                                                                                                                                                                                                                                                                                                     | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C178  | —     | Pearl Pink                                                                                                                                                                                                                                                                                                                                                      | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C182  | —     | Silver                                                                                                                                                                                                                                                                                                                                                          | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C125  | —     | White                                                                                                                                                                                                                                                                                                                                                           | NC                                                                 |            |
| <input type="checkbox"/>                                                             | 400C170  | —     | American Flag                                                                                                                                                                                                                                                                                                                                                   | \$602.00                                                           |            |

| BASE                                                                                                                                                                                                                     |       |                                                         |          |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|---------------------------------------------------------|----------|
| OPTION                                                                                                                                                                                                                   | HCPC  | DESCRIPTION                                             | RETAIL   |
| <b>Color Accents</b> <i>Note: Color Accents and Shroud Colors may not match.</i>                                                                                                                                         |       |                                                         | \$       |
| <input type="checkbox"/> 400EZ20                                                                                                                                                                                         | —     | Apple Green                                             | NC       |
| <input type="checkbox"/> 400EZ22                                                                                                                                                                                         | —     | Black                                                   | NC       |
| <input type="checkbox"/> 400EZ21                                                                                                                                                                                         | —     | Blue                                                    | NC       |
| <input type="checkbox"/> 400EZ23                                                                                                                                                                                         | —     | Orange                                                  | NC       |
| <input type="checkbox"/> 400EZ24                                                                                                                                                                                         | —     | Red                                                     | NC       |
| <input type="checkbox"/> 400EZ19                                                                                                                                                                                         | —     | White                                                   | NC       |
| <input type="checkbox"/> 400EZ25                                                                                                                                                                                         | —     | Silver                                                  | NC       |
| <b>Service Decals</b>                                                                                                                                                                                                    |       |                                                         | \$       |
| <input type="checkbox"/> 400D10                                                                                                                                                                                          | —     | Army                                                    | NC       |
| <input type="checkbox"/> 400D11                                                                                                                                                                                          | —     | Air Force                                               | NC       |
| <input type="checkbox"/> 400D13                                                                                                                                                                                          | —     | Marines                                                 | NC       |
| <input type="checkbox"/> 400D14                                                                                                                                                                                          | —     | Navy                                                    | NC       |
| <input type="checkbox"/> 400D15                                                                                                                                                                                          | —     | Purple Heart                                            | NC       |
| <b>ISO Transit Tie Downs</b>                                                                                                                                                                                             |       |                                                         | \$       |
| <input type="checkbox"/> 400T050                                                                                                                                                                                         | —     | ISO Transit Package                                     | NC       |
| <b>LED Lights</b>                                                                                                                                                                                                        |       |                                                         | \$       |
| <input type="checkbox"/> 400LI08                                                                                                                                                                                         | —     | LED Fender Lights                                       | \$566.00 |
| <input type="checkbox"/> 400LI09                                                                                                                                                                                         | —     | LED Fender Lights <i>(Available w/Multi-Power only)</i> | NC       |
| <input type="checkbox"/> 400LI02                                                                                                                                                                                         | —     | Omit LED Lights                                         | NC       |
| <b>Seat Heights</b>                                                                                                                                                                                                      |       |                                                         | \$       |
| <b>Seat Height</b> <i>Measured from front of seat pan to floor with 0° frame angle. Not all seat heights are available with all seat depths and power seating. Please see online configurator for full availability.</i> |       |                                                         |          |
| <input type="checkbox"/> 400HT103                                                                                                                                                                                        | —     | 16.5"                                                   | NC       |
| <input type="checkbox"/> 400HT55                                                                                                                                                                                         | —     | 17"                                                     | NC       |
| <input type="checkbox"/> 400HT104                                                                                                                                                                                        | —     | 17.5"                                                   | NC       |
| <input type="checkbox"/> 400HT56                                                                                                                                                                                         | —     | 18"                                                     | NC       |
| <input type="checkbox"/> 400HT105                                                                                                                                                                                        | —     | 18.5"                                                   | NC       |
| <input type="checkbox"/> 400HT57                                                                                                                                                                                         | —     | 19"                                                     | NC       |
| <b>Batteries and Chargers</b>                                                                                                                                                                                            |       |                                                         | \$       |
| <b>Group 22 Battery</b>                                                                                                                                                                                                  |       |                                                         |          |
| <input type="checkbox"/> 400AG35                                                                                                                                                                                         | E2361 | GP22 / 50Ah Flat Top Pair                               | \$680.00 |
| <b>Omit Batteries</b>                                                                                                                                                                                                    |       |                                                         |          |
| <input type="checkbox"/> 400AG24                                                                                                                                                                                         | —     | Battery Harnesses Only <i>(Omit Batteries)</i>          | NC       |
| <b>Chargers</b>                                                                                                                                                                                                          |       |                                                         |          |
| <input type="checkbox"/> 400BJ30                                                                                                                                                                                         | —     | 10A Dual Mode (Compact)                                 | \$110.00 |
| <input type="checkbox"/> 400BJ19                                                                                                                                                                                         | —     | 8A Dual Mode (Compact)                                  | NC       |

| POWERED SEATING                                                                            |       |                                                                                                                          |             |
|--------------------------------------------------------------------------------------------|-------|--------------------------------------------------------------------------------------------------------------------------|-------------|
| OPTION                                                                                     | HCPC  | DESCRIPTION                                                                                                              | RETAIL      |
| Power Seating Options                                                                      |       |                                                                                                                          | \$          |
| <b>Seat Elevator</b>                                                                       |       |                                                                                                                          |             |
| <input type="checkbox"/> 400PI07                                                           | E2298 | Seat Elevator, Powered, 300 mm / 12"                                                                                     | \$3,791.00  |
| <input type="checkbox"/> 400PI08                                                           | E2298 | C-Me® Seat Elevator, Powered, 300 mm / 12" <i>(C-Me speed up to 3 mph) (N/A with VR-2)</i>                               | \$3,791.00  |
| <b>Seat Tilt</b>                                                                           |       |                                                                                                                          |             |
| <input type="checkbox"/> 400PA48                                                           | E1002 | Posterior 50° Tilt                                                                                                       | \$8,458.00  |
| <b>Seat Recline</b>                                                                        |       |                                                                                                                          |             |
| <input type="checkbox"/> 400PA37                                                           | E1004 | 150° Powered Recline w/ 4" Shear Reduction                                                                               | \$8,458.00  |
| <b>Seat Tilt/Recline</b>                                                                   |       |                                                                                                                          |             |
| <input type="checkbox"/> 400PA41                                                           | E1007 | Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction                                                       | \$13,709.00 |
| <b>Power Recline Backrest Preset Home Angle</b>                                            |       |                                                                                                                          |             |
| <input type="checkbox"/> 400QK1                                                            | —     | 90°                                                                                                                      | NC          |
| <input type="checkbox"/> 400QK2                                                            | —     | 95° (Default Position)                                                                                                   | NC          |
| Pre Set Seat Frame Angle <i>(No Power Seating Options, Increases Seat to Floor Height)</i> |       |                                                                                                                          | \$          |
| <input type="checkbox"/> 400FA10                                                           | —     | 0° <i>(Adds 0" to STFH)</i>                                                                                              | NC          |
| <input type="checkbox"/> 400FA3                                                            | —     | 3° <i>(Adds 1" to STFH)</i>                                                                                              | NC          |
| Backrest Frame Type                                                                        |       |                                                                                                                          | \$          |
| <input type="checkbox"/> 400BF16                                                           | —     | Tool Angle Adjustable Back, 90°-120°                                                                                     | NC          |
| <input type="checkbox"/> 400BF14                                                           | E1225 | User Adjustable Manual Semi-Recline Back, 90°-130°, Right <i>(Adjustment Lever Mounted on Opposite Side of Joystick)</i> | \$682.00    |
| <input type="checkbox"/> 400BF15                                                           | E1225 | User Adjustable Manual Semi-Recline Back, 90°-130°, Left <i>(Adjustment Lever Mounted on Opposite Side of Joystick)</i>  | \$682.00    |
| Backrest Frame Option                                                                      |       |                                                                                                                          | \$          |
| <input type="checkbox"/> 400BF13                                                           | —     | Flip Down, Quick Release                                                                                                 | \$147.00    |
| Back Cushion Type                                                                          |       |                                                                                                                          | \$          |
| <input type="checkbox"/> 400BT134                                                          | —     | SEDEO® Back Cushion <i>(Choose Seating in Section A)</i>                                                                 | NC          |
| <input type="checkbox"/> 400BT138                                                          | —     | JAY® Syncra™ <i>(Choose Seating in Section B)</i>                                                                        | NC          |
| <input type="checkbox"/> 400BT135                                                          | —     | JAY® J3® <i>(Choose Seating in Section C)</i>                                                                            | NC          |
| <input type="checkbox"/> 400BT137                                                          | —     | SEDEO® PRO™ Aftermarket Back Adapter <i>(Height Adjustable Round Tubes for Traditional Mounting)</i>                     | \$525.00    |
| <input type="checkbox"/> 400BT133                                                          | —     | OMIT <i>(For personalized Mounting. No Headrest Mount)</i>                                                               | NC          |

## POWERED SEATING

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

## A - SEDEO® Backrest Dimensions

\$

## Back Width

|                          |         |   |             |    |
|--------------------------|---------|---|-------------|----|
| <input type="checkbox"/> | 400BN28 | — | 400mm / 16" | NC |
| <input type="checkbox"/> | 400BN29 | — | 460mm / 18" | NC |
| <input type="checkbox"/> | 400BN30 | — | 500mm / 20" | NC |
| <input type="checkbox"/> | 400BN31 | — | 560mm / 22" | NC |

## Back Shell Height

## Tool Angle Adjust.

## 150° Recline

|                          |         |   |             |         |         |                       |    |
|--------------------------|---------|---|-------------|---------|---------|-----------------------|----|
| <input type="checkbox"/> | 400B255 | — | 460mm / 18" | 22"-23" | 23"-25" | Installed Back Height | NC |
| <input type="checkbox"/> | 400B256 | — | 560mm / 22" | 24"-27" | 26"-29" |                       | NC |
| <input type="checkbox"/> | 400B257 | — | 660mm / 26" | 28"-30" | 30"     |                       | NC |

## Installed Back Height - (From Pan to Top of Back)

|                          |         |   |     |    |
|--------------------------|---------|---|-----|----|
| <input type="checkbox"/> | 400JR32 | — | 22" | NC |
| <input type="checkbox"/> | 400JR33 | — | 23" | NC |
| <input type="checkbox"/> | 400JR34 | — | 24" | NC |
| <input type="checkbox"/> | 400JR35 | — | 25" | NC |
| <input type="checkbox"/> | 400JR36 | — | 26" | NC |
| <input type="checkbox"/> | 400JR37 | — | 27" | NC |
| <input type="checkbox"/> | 400JR38 | — | 28" | NC |
| <input type="checkbox"/> | 400JR39 | — | 29" | NC |
| <input type="checkbox"/> | 400JR40 | — | 30" | NC |

## B - JAY® Syncra™ Seating

B - JAY® Syncra™ Back  P. 21

\$

|                          |         |       |                                                     |            |
|--------------------------|---------|-------|-----------------------------------------------------|------------|
| <input type="checkbox"/> | 400JP14 | E2620 | JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth | \$1,138.00 |
| <input type="checkbox"/> | 400JP15 | E2621 | JAY® Syncra™ Back 22" Wide with 3" Contour Depth    | \$1,255.00 |

## B - JAY® Syncra™ Back Width

\$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

|                          |         |   |                                             |    |
|--------------------------|---------|---|---------------------------------------------|----|
| <input type="checkbox"/> | 400JQ07 | — | 16" Wide (Adjustment range 16 to 17" Wide)* | NC |
| <input type="checkbox"/> | 400JQ08 | — | 18" Wide (Adjustment range 18 to 19" Wide)* | NC |
| <input type="checkbox"/> | 400JQ09 | — | 20" Wide (Adjustment range 20 to 21" Wide)* | NC |
| <input type="checkbox"/> | 400JQ10 | — | 22" Wide (Adjustment range 22 to 23" Wide)  | NC |

\*Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.

**POWERED SEATING**

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

**B - JAY® Syncra™ Back Shell Options and Height Ranges**

\$

JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown based on the back frame type.

**JAY® Syncra™ Shell Height**

|                          |         |   |                                         |    |
|--------------------------|---------|---|-----------------------------------------|----|
| <input type="checkbox"/> | 400JR42 | – | 16" Tall Shell (22" to 24" Back Height) | NC |
| <input type="checkbox"/> | 400JR43 | – | 18" Tall Shell (24" to 26" Back Height) | NC |
| <input type="checkbox"/> | 400JR44 | – | 20" Tall Shell (26" to 28" Back Height) | NC |
| <input type="checkbox"/> | 400JR45 | – | 22" Tall Shell (28" to 30" Back Height) | NC |

**JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)**

|                          |         |   |                        |    |
|--------------------------|---------|---|------------------------|----|
| <input type="checkbox"/> | 400JR50 | – | 22" (16" Shell Height) | NC |
| <input type="checkbox"/> | 400JR51 | – | 23" (16" Shell Height) | NC |
| <input type="checkbox"/> | 400JR52 | – | 24" (16" Shell Height) | NC |
| <input type="checkbox"/> | 400JR53 | – | 24" (18" Shell Height) | NC |
| <input type="checkbox"/> | 400JR54 | – | 25" (18" Shell Height) | NC |
| <input type="checkbox"/> | 400JR55 | – | 26" (18" Shell Height) | NC |
| <input type="checkbox"/> | 400JR56 | – | 26" (20" Shell Height) | NC |
| <input type="checkbox"/> | 400JR57 | – | 27" (20" Shell Height) | NC |
| <input type="checkbox"/> | 400JR58 | – | 28" (20" Shell Height) | NC |
| <input type="checkbox"/> | 400JR67 | – | 28" (22" Shell Height) | NC |
| <input type="checkbox"/> | 400JR68 | – | 29" (22" Shell Height) | NC |
| <input type="checkbox"/> | 400JR69 | – | 30" (22" Shell Height) | NC |

**B - JAY® Syncra™ Back Upholstery**

\$

|                          |         |   |                             |    |
|--------------------------|---------|---|-----------------------------|----|
| <input type="checkbox"/> | 400BU96 | – | Syncra™ Premium Leatherette | NC |
| <input type="checkbox"/> | 400BU97 | – | Syncra™ Stretch             | NC |

**C - JAY® J3® Back****C - JAY® J3® Back - Contour**

\$

|                          |         |       |                                                              |          |
|--------------------------|---------|-------|--------------------------------------------------------------|----------|
| <input type="checkbox"/> | 400JP01 | E2613 | J3® PA Back for 16"-20" Wide                                 | \$625.00 |
| <input type="checkbox"/> | 400JP02 | E2614 | J3® PA Back for 22" Wide                                     | \$928.00 |
| <input type="checkbox"/> | 400JP11 | E2615 | J3® PL Back for 16"-20" Wide                                 | \$729.00 |
| <input type="checkbox"/> | 400JP13 | E2620 | J3® PD Back for 16"-20" Wide (18" and 21" Shell Height Only) | \$875.00 |

**C - JAY® J3® Back Width**

\$

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

|                          |         |   |     |    |
|--------------------------|---------|---|-----|----|
| <input type="checkbox"/> | 400JQ02 | – | 16" | NC |
| <input type="checkbox"/> | 400JQ03 | – | 18" | NC |
| <input type="checkbox"/> | 400JQ04 | – | 20" | NC |
| <input type="checkbox"/> | 400JQ05 | – | 22" | NC |

## POWERED SEATING

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

## C - JAY® J3® Back Shell Options and Height Ranges

\$

JAY® J3® shells are mounted with a minimum 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.

| JAY® J3® Shell Height    |         |   | Tool Angle Adjust.                               | 150° Recline | Installed Back Height |    |
|--------------------------|---------|---|--------------------------------------------------|--------------|-----------------------|----|
| <input type="checkbox"/> | 400JR41 | — | 15" tall shell (Upper Thoracic Med) for 22" Wide | 22"-23"      |                       | NC |
| <input type="checkbox"/> | 400JR14 | — | 16.5" tall shell (Upper Thoracic Short)          | 21"-23"      |                       | NC |
| <input type="checkbox"/> | 400JR15 | — | 18" tall shell (Upper Thoracic Med)              | 22"-24"      |                       | NC |
| <input type="checkbox"/> | 400JR16 | — | 19.5" tall shell (Upper Thoracic Tall)           | 24"-26"      |                       | NC |
| <input type="checkbox"/> | 400JR18 | — | 18" tall shell (Shoulder Height) for 22" wide    | 24"          |                       | NC |
| <input type="checkbox"/> | 400JR17 | — | 21" tall shell (Shoulder Height)                 | 24"-26"      |                       | NC |
| <input type="checkbox"/> | 400JR27 | — | 22.5" tall shell (Shoulder Height)               | 26"-28"      |                       | NC |
| <input type="checkbox"/> | 400JR28 | — | 24" tall shell (Shoulder Height)                 | 27"-29"      |                       | NC |

## JAY® J3® Installed Back Height (From Pan to Top of Back)

|                          |         |   |     |  |    |
|--------------------------|---------|---|-----|--|----|
| <input type="checkbox"/> | 400JR29 | — | 21" |  | NC |
| <input type="checkbox"/> | 400JR19 | — | 22" |  | NC |
| <input type="checkbox"/> | 400JR20 | — | 23" |  | NC |
| <input type="checkbox"/> | 400JR21 | — | 24" |  | NC |
| <input type="checkbox"/> | 400JR22 | — | 25" |  | NC |
| <input type="checkbox"/> | 400JR23 | — | 26" |  | NC |
| <input type="checkbox"/> | 400JR24 | — | 27" |  | NC |
| <input type="checkbox"/> | 400JR25 | — | 28" |  | NC |
| <input type="checkbox"/> | 400JR30 | — | 29" |  | NC |
| <input type="checkbox"/> | 400JR31 | — | 30" |  | NC |

## C - JAY® Anterior Thoracic Support

\$

Classic-Structure Pull Down Harness (4 point strap, Only available with JAY® backrest)

|                          |        |       |                  |          |
|--------------------------|--------|-------|------------------|----------|
| <input type="checkbox"/> | 400K10 | E0960 | Medium (ATS8302) | \$177.00 |
| <input type="checkbox"/> | 400K11 | E0960 | Large (ATS8303)  | \$177.00 |

## C - JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)

\$

## JAY® Lateral Thoracic Support Pads

|                          |         |       |                                                                |          |
|--------------------------|---------|-------|----------------------------------------------------------------|----------|
| <input type="checkbox"/> | 400JW35 | E0956 | Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right | \$147.00 |
| <input type="checkbox"/> | 400JW38 | E0956 | Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left  | \$147.00 |
| <input type="checkbox"/> | 400JW36 | E0956 | Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right | \$147.00 |
| <input type="checkbox"/> | 400JW39 | E0956 | Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left   | \$147.00 |
| <input type="checkbox"/> | 400JW37 | E0956 | Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right | \$147.00 |
| <input type="checkbox"/> | 400JW40 | E0956 | Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left  | \$147.00 |

**POWERED SEATING**

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

**JAY® Lateral Thoracic Support Pad Upholstery**

|                          |         |   |                             |    |
|--------------------------|---------|---|-----------------------------|----|
| <input type="checkbox"/> | 400NU01 | – | Syncra™ Premium Leatherette | NC |
| <input type="checkbox"/> | 400NU02 | – | Syncra™ Stretch             | NC |

**JAY® Lateral Thoracic Swing-Away Mount**

|                          |         |                |                                             |          |
|--------------------------|---------|----------------|---------------------------------------------|----------|
| <input type="checkbox"/> | 400JW41 | E1028<br>E1034 | Lateral Thoracic Swing Away Bracket - Right | \$292.00 |
| <input type="checkbox"/> | 400JW42 | E1028<br>E1034 | Lateral Thoracic Swing Away Bracket - Left  | \$292.00 |

**C - JAY® Syncra™ Fixed Lateral Support Accessory**

\$

**JAY® Fixed Lateral Support Accessory Hardware** (JAY® Fixed Lateral Support is height adjustable)

|                          |         |       |                                                                            |          |
|--------------------------|---------|-------|----------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 400JW33 | E0956 | Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right | \$147.00 |
| <input type="checkbox"/> | 400JW34 | E0956 | Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left  | \$147.00 |

**Headrests (Requires Pad, Upholstery, & Mount)**

\$

**Headrest Pad**

|                          |         |       |                                                                                                                                                                  |          |
|--------------------------|---------|-------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 400HE36 | E0955 | SEDEO® Pad Small with Cover (LxH) 240 x 150 mm / 10" x 6"<br><i>Can only be used with Multi Position Headrest Mount with quick-release and SEDEO® upholstery</i> | \$188.00 |
| <input type="checkbox"/> | 400HE38 | E0955 | Whitmyer® Plush Headrest Pad 203 mm / 8"                                                                                                                         | \$239.00 |
| <input type="checkbox"/> | 400HE39 | E0955 | Whitmyer® Plush Headrest Pad 254 mm / 10"                                                                                                                        | \$239.00 |
| <input type="checkbox"/> | 400HE44 | E0955 | Whitmyer® Plush Headrest Pad 356 mm / 14"                                                                                                                        | \$250.00 |
| <input type="checkbox"/> | 400HE46 | E0955 | Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable 7" to 16" Wide x 4" Tall)                                                                               | \$367.00 |
| <input type="checkbox"/> | 400HE48 | –     | OMIT Headrest Pad                                                                                                                                                | NC       |

**Headrest Pad Upholstery**

|                          |         |   |                         |    |
|--------------------------|---------|---|-------------------------|----|
| <input type="checkbox"/> | 400NU10 | – | SEDEO® Upholstery       | NC |
| <input type="checkbox"/> | 400NU11 | – | Whitmyer® Lycra/Stretch | NC |

**Headrest Mount**

|                          |         |                |                                                                                          |          |
|--------------------------|---------|----------------|------------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 400MH32 | E1028<br>E1033 | SEDEO® Multi-Position Headrest Mount with Quick Release (Only available with SEDEO® Pad) | \$234.00 |
| <input type="checkbox"/> | 400MH33 | E1028<br>E1033 | Whitmyer® Linx Headrest Mount with Quick Release                                         | \$192.00 |
| <input type="checkbox"/> | 400MH34 | E1028<br>E1033 | Whitmyer® Linx-2 Headrest Mount with Quick Release                                       | \$234.00 |
| <input type="checkbox"/> | 400MH35 | E1028<br>E1033 | Whitmyer® AXYS Headrest Mount with Quick Release                                         | \$123.00 |
| <input type="checkbox"/> | 400MH17 | –              | OMIT Headrest Mount                                                                      | NC       |

**Push Handles**

\$

|                          |         |   |                                |          |
|--------------------------|---------|---|--------------------------------|----------|
| <input type="checkbox"/> | 400PH27 | – | Height Adjustable Push Handles | \$350.00 |
|--------------------------|---------|---|--------------------------------|----------|

**Seat Frame Rails**

\$

|                          |         |   |                                                     |    |
|--------------------------|---------|---|-----------------------------------------------------|----|
| <input type="checkbox"/> | 400SN14 | – | Seat Rail Short, for Effective Seat Depth 14" - 18" | NC |
| <input type="checkbox"/> | 400SN15 | – | Seat Rail Long, for Effective Seat Depth 17" - 22"  | NC |

## POWERED SEATING

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

Seat Cushion Type \$

|                          |         |   |                              |    |
|--------------------------|---------|---|------------------------------|----|
| <input type="checkbox"/> | 400ST81 | – | JAY® Syncra™ (See Section D) | NC |
| <input type="checkbox"/> | 400ST80 | – | JAY® (See Section E)         | NC |
| <input type="checkbox"/> | 400ST75 | – | Omit                         | NC |

Seat Size \$

## Seat Pan Width

|                          |        |   |                                                       |    |
|--------------------------|--------|---|-------------------------------------------------------|----|
| <input type="checkbox"/> | 400R53 | – | 16"                                                   | NC |
| <input type="checkbox"/> | 400R54 | – | 17"                                                   | NC |
| <input type="checkbox"/> | 400R55 | – | 18"                                                   | NC |
| <input type="checkbox"/> | 400R56 | – | 19"                                                   | NC |
| <input type="checkbox"/> | 400R57 | – | 20"                                                   | NC |
| <input type="checkbox"/> | 400R58 | – | 21" <i>(Not available with JAY® Syncra™ Cushions)</i> | NC |
| <input type="checkbox"/> | 400R59 | – | 22"                                                   | NC |

Seat Pan Depth *(Front of seat pan to front of back cushion – effective depth)*

|                          |         |   |                                                       |    |
|--------------------------|---------|---|-------------------------------------------------------|----|
| <input type="checkbox"/> | 400RF52 | – | 14"                                                   | NC |
| <input type="checkbox"/> | 400RF53 | – | 15"                                                   | NC |
| <input type="checkbox"/> | 400RF51 | – | 16"                                                   | NC |
| <input type="checkbox"/> | 400RF45 | – | 17"                                                   | NC |
| <input type="checkbox"/> | 400RF46 | – | 18"                                                   | NC |
| <input type="checkbox"/> | 400RF47 | – | 19"                                                   | NC |
| <input type="checkbox"/> | 400RF48 | – | 20"                                                   | NC |
| <input type="checkbox"/> | 400RF49 | – | 21" <i>(Not available with JAY® Syncra™ Cushions)</i> | NC |
| <input type="checkbox"/> | 400RF50 | – | 22"                                                   | NC |

D - JAY® Syncra™ Cushion  P.21 (JAY Your Way modifications available through Customer Service (800) 333-4000) \$

|                          |         |       |                                                                                               |          |
|--------------------------|---------|-------|-----------------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 400CU89 | E2622 | JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology <i>(16" to 20" wide)</i> | \$630.00 |
| <input type="checkbox"/> | 400CU90 | E2623 | JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology <i>(22" wide)</i>        | \$729.00 |

## JAY® Syncra™ Cushion Cover Upholstery

|                          |         |   |                             |    |
|--------------------------|---------|---|-----------------------------|----|
| <input type="checkbox"/> | 400YR73 | – | Syncra™ Premium Leatherette | NC |
| <input type="checkbox"/> | 400YR74 | – | Syncra™ Stretch             | NC |

E - JAY® Seat Cushion \$

|                          |         |       |                                             |          |
|--------------------------|---------|-------|---------------------------------------------|----------|
| <input type="checkbox"/> | 400CU59 | E2607 | JAY® Ion® <i>(Sizes 16" to 21" Wide)</i>    | \$460.00 |
| <input type="checkbox"/> | 400CU60 | E2608 | JAY® Ion® <i>(Size 22" Wide)</i>            | \$578.00 |
| <input type="checkbox"/> | 400CU61 | E2607 | JAY® Union® <i>(Sizes 16" to 21" Wide)</i>  | \$520.00 |
| <input type="checkbox"/> | 400CU62 | E2608 | JAY® Union® <i>(Size 22" Wide)</i>          | \$635.00 |
| <input type="checkbox"/> | 400CU64 | E2622 | JAY® Fusion® <i>(Sizes 16" to 21" Wide)</i> | \$514.00 |
| <input type="checkbox"/> | 400CU65 | E2623 | JAY® Fusion® <i>(Size 22" Wide)</i>         | \$613.00 |

## POSITIONING SUPPORTS &amp; ACCESSORIES

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

## Knee and Thigh Swing-Away Support Pads

\$

Lateral Knee Support Pads and Hardware *(must pick pad and mount)*

|                          |         |       |                                            |          |
|--------------------------|---------|-------|--------------------------------------------|----------|
| <input type="checkbox"/> | 400XY05 | E0953 | Lateral Knee Support Pad, Right (4" round) | \$205.00 |
| <input type="checkbox"/> | 400XY06 | E0953 | Lateral Knee Support Pad, Left (4" round)  | \$205.00 |
| <input type="checkbox"/> | 400XX08 | E1028 | Lateral Swing-Away Bracket Knee - Right    | \$268.00 |
| <input type="checkbox"/> | 400XX09 | E1028 | Lateral Swing-Away Bracket Knee - Left     | \$268.00 |

Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*

|                          |         |       |                                                      |          |
|--------------------------|---------|-------|------------------------------------------------------|----------|
| <input type="checkbox"/> | 400KF22 | E0953 | Lateral Thigh Support Pad Short, Right (5" H x 5" L) | \$147.00 |
| <input type="checkbox"/> | 400KF23 | E0953 | Lateral Thigh Support Pad Short, Left (5" H x 5" L)  | \$147.00 |
| <input type="checkbox"/> | 400KF24 | E0953 | Lateral Thigh Support Pad Long, Right (5" H x 7" L)  | \$147.00 |
| <input type="checkbox"/> | 400KF25 | E0953 | Lateral Thigh Support Pad Long, Left (5" H x 7" L)   | \$147.00 |
| <input type="checkbox"/> | 400KF26 | E1028 | Lateral Thigh Removable Mount - Right                | \$292.00 |
| <input type="checkbox"/> | 400KF27 | E1028 | Lateral Thigh Removable Mount - Left                 | \$292.00 |

JAY® Lateral Thigh Supports *(Must pick pad, upholstery, and mount)*

\$

## JAY® Lateral Thigh Support Pads

|                          |         |       |                                                           |          |
|--------------------------|---------|-------|-----------------------------------------------------------|----------|
| <input type="checkbox"/> | 400KF16 | E0953 | Lateral Thigh Support Pad Standard, Right (8.25" x 3.25") | \$147.00 |
| <input type="checkbox"/> | 400KF17 | E0953 | Lateral Thigh Support Pad Standard, Left (8.25" x 3.25")  | \$147.00 |
| <input type="checkbox"/> | 400KF18 | E0953 | Lateral Thigh Support Pad Large, Right (12" x 3.25")      | \$147.00 |
| <input type="checkbox"/> | 400KF19 | E0953 | Lateral Thigh Support Pad Large, Left (12" x 3.25")       | \$147.00 |

## JAY® Lateral Thigh Support Pad Upholstery

|                          |         |   |                             |    |
|--------------------------|---------|---|-----------------------------|----|
| <input type="checkbox"/> | 400NU07 | — | Syncra™ Premium Leatherette | NC |
| <input type="checkbox"/> | 400NU08 | — | Syncra™ Stretch             | NC |

## JAY® Lateral Thigh Removable Mount

|                          |         |       |                                       |          |
|--------------------------|---------|-------|---------------------------------------|----------|
| <input type="checkbox"/> | 400KF20 | E1028 | Lateral Thigh Removable Mount - Right | \$292.00 |
| <input type="checkbox"/> | 400KF21 | E1028 | Lateral Thigh Removable Mount - Left  | \$292.00 |

## Armrest Structure

\$

|                          |         |   |                                                                                                               |          |
|--------------------------|---------|---|---------------------------------------------------------------------------------------------------------------|----------|
| <input type="checkbox"/> | 400EG41 | — | Reclining, Flip Back, Short <i>(Best used with shorter seat depths)</i> (Height Adjustable 9" to 15")         | NC       |
| <input type="checkbox"/> | 400EG36 | — | Reclining, Flip Back with 2-point post locking arm (Height Adjustable 11" to 15" x 0.625")                    | \$292.00 |
| <input type="checkbox"/> | 400EG52 | — | Reclining, Flip Back, Dual Link Armrest (Height Adjustable 8" to 15" including length and lateral adjustment) | NC       |
| <input type="checkbox"/> | 400EG32 | — | Single Post, Flip Back <i>(N/A with thigh support)</i> (Height Adjustable 10" to 15")                         | NC       |
| <input type="checkbox"/> | 400EG50 | — | Dual Post, Flip Back, Removable (Height Adjustable)                                                           | \$225.00 |

## POSITIONING SUPPORTS &amp; ACCESSORIES

| OPTION | HCPC | DESCRIPTION | RETAIL |
|--------|------|-------------|--------|
|--------|------|-------------|--------|

## Arm Pads

\$

## SEDEO® Armrest Pad

|                          |         |   |                                    |    |
|--------------------------|---------|---|------------------------------------|----|
| <input type="checkbox"/> | 400AB43 | — | 300 mm / 12", Integral Foam, Right | NC |
| <input type="checkbox"/> | 400AE43 | — | 300 mm / 12", Integral Foam, Left  | NC |
| <input type="checkbox"/> | 400AB44 | — | 400 mm / 16", Integral Foam, Right | NC |
| <input type="checkbox"/> | 400AE44 | — | 400 mm / 16", Integral Foam, Left  | NC |

## Dual Post Armrest Pad

|                          |         |   |                                        |    |
|--------------------------|---------|---|----------------------------------------|----|
| <input type="checkbox"/> | 400AB74 | — | 250 mm / 10", Classic Desk Foam, Right | NC |
| <input type="checkbox"/> | 400AE74 | — | 250 mm / 10", Classic Desk Foam, Left  | NC |
| <input type="checkbox"/> | 400AB76 | — | 250 mm / 10", Waterfall Desk, Right    | NC |
| <input type="checkbox"/> | 400AE76 | — | 250 mm / 10", Waterfall Desk, Left     | NC |
| <input type="checkbox"/> | 400AB75 | — | 350 mm / 14", Classic Full Foam, Right | NC |
| <input type="checkbox"/> | 400AE75 | — | 350 mm / 14", Classic Full Foam, Left  | NC |
| <input type="checkbox"/> | 400AB77 | — | 350 mm / 14", Waterfall Full, Right    | NC |
| <input type="checkbox"/> | 400AE77 | — | 350 mm / 14", Waterfall Full, Left     | NC |

## Comfort Arm Pad (Foam)

|                          |         |   |               |         |
|--------------------------|---------|---|---------------|---------|
| <input type="checkbox"/> | 400AB46 | — | 2" Wide Right | \$88.00 |
| <input type="checkbox"/> | 400AE46 | — | 2" Wide Left  | \$88.00 |

## Comfort Arm Pad (Gel)

|                          |         |   |               |          |
|--------------------------|---------|---|---------------|----------|
| <input type="checkbox"/> | 400AB48 | — | 2" Wide Right | \$175.00 |
| <input type="checkbox"/> | 400AE48 | — | 2" Wide Left  | \$175.00 |

## Ottobock® Arm Troughs

|                          |         |   |                                          |          |
|--------------------------|---------|---|------------------------------------------|----------|
| <input type="checkbox"/> | 400AE06 | — | Ottobock Arm Trough, Med (12"), Left     | \$136.00 |
| <input type="checkbox"/> | 400AB06 | — | Ottobock Arm Trough, Med (12"), Right    | \$136.00 |
| <input type="checkbox"/> | 400AE07 | — | Ottobock Arm Trough, Large (14"), Left   | \$136.00 |
| <input type="checkbox"/> | 400AB07 | — | Ottobock Arm Trough, Large (14"), Right  | \$136.00 |
| <input type="checkbox"/> | 400OL01 | — | Ottobock Hand Pad, Med (7"), Left        | \$59.00  |
| <input type="checkbox"/> | 400OR01 | — | Ottobock Hand Pad, Med (7"), Right       | \$59.00  |
| <input type="checkbox"/> | 400OL02 | — | Ottobock Hand Pad, Large (9"), Left      | \$59.00  |
| <input type="checkbox"/> | 400OR02 | — | Ottobock Hand Pad, Large (9"), Right     | \$59.00  |
| <input type="checkbox"/> | 400AE78 | — | Ottobock 1 Piece Arm Trough (18"), Left  | \$190.00 |
| <input type="checkbox"/> | 400AB78 | — | Ottobock 1 Piece Arm Trough (18"), Right | \$190.00 |

## POSITIONING SUPPORTS & ACCESSORIES

| OPTION                                                                                                              | HCPC  | DESCRIPTION                                                                                                             | RETAIL     |
|---------------------------------------------------------------------------------------------------------------------|-------|-------------------------------------------------------------------------------------------------------------------------|------------|
| <b>Elbow Support</b>                                                                                                |       |                                                                                                                         | <b>\$</b>  |
| <input type="checkbox"/> 4003A02                                                                                    | —     | Elbow Support Pad, RH 3" x 10" <i>(Includes side support)</i>                                                           | \$179.00   |
| <input type="checkbox"/> 4003A03                                                                                    | —     | Elbow Support Pad, LH 3" x 10" <i>(Includes side support)</i>                                                           | \$179.00   |
| <input type="checkbox"/> 4003A06                                                                                    | —     | Elbow Support Pad, RH 3" x 3"                                                                                           | \$149.00   |
| <input type="checkbox"/> 4003A07                                                                                    | —     | Elbow Support Pad, LH 3" x 3"                                                                                           | \$149.00   |
| <input type="checkbox"/> 4003A04                                                                                    | E1028 | Swing-Away Elbow Support Hardware, RH                                                                                   | \$245.00   |
| <input type="checkbox"/> 4003A05                                                                                    | E1028 | Swing-Away Elbow Support Hardware, LH                                                                                   | \$245.00   |
| <b>Side Guards (Pair)</b>                                                                                           |       |                                                                                                                         | <b>\$</b>  |
| <input type="checkbox"/> 400SG44                                                                                    | E0953 | Side Guards Integral Foam (11"L x 5.5"H) <i>(N/A with Reclining Armrests)</i>                                           | \$239.00   |
| <b>Fixed Centermount Legrest</b>                                                                                    |       |                                                                                                                         | <b>\$</b>  |
| <input type="checkbox"/> 400H184                                                                                    | —     | Fixed Centermount Legrest, Short <i>(Includes Calf Pads) (11.5" to 14.5" Leg Length)</i>                                | NC         |
| <input type="checkbox"/> 400H185                                                                                    | —     | Fixed Centermount Legrest, Long <i>(Includes Calf Pads) (13.5" to 18.5" Leg Length)</i>                                 | NC         |
| <input type="checkbox"/> 400H136                                                                                    | —     | Angle Adjustable Centermount Legrest <i>(11.5" to 18.5" Leg Length)</i>                                                 | \$525.00   |
| <b>Power Centermount Legrest</b>                                                                                    |       |                                                                                                                         | <b>\$</b>  |
| <input type="checkbox"/> 400H109                                                                                    | E1012 | SEDEO® Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension <i>(Powered Articulation) (3" to 17.5" Leg Length)</i> | \$3,791.00 |
| <b>Centermount Legrest Accessories</b>                                                                              |       |                                                                                                                         | <b>\$</b>  |
| <b>Power Centermount Footplates</b>                                                                                 |       |                                                                                                                         |            |
| <input type="checkbox"/> 400F127                                                                                    | K0040 | Dual, Height, Width and Angle Adjustable, Standard <i>(14"W x 10"D)</i>                                                 | \$197.00   |
| <input type="checkbox"/> 400F128                                                                                    | K0040 | Dual, Height, Width and Angle Adjustable, Large <i>(16"W x 12"D)</i>                                                    | \$197.00   |
| <b>Footplate Accessories (Available with Power Centermounts)</b>                                                    |       |                                                                                                                         |            |
| <input type="checkbox"/> 400DI68                                                                                    | K0108 | Extra Large Centermount Footplate <i>(Single Footplate Overlay 16"W x 12"D)</i>                                         | \$199.00   |
| <input type="checkbox"/> 400DI69                                                                                    | K0108 | Footplate Connector Link <i>(Attaches footplates together)</i>                                                          | \$97.00    |
| <input type="checkbox"/> 400DI57                                                                                    | K0108 | Gel Footpads <i>(Pair, LH and RH)</i>                                                                                   | \$403.00   |
| <b>Angle Adjustable / Fixed Centermount Footplates</b>                                                              |       |                                                                                                                         |            |
| <input type="checkbox"/> 400F45                                                                                     | —     | Single/Dual Footplate                                                                                                   | NC         |
| <b>Power Centermount Calf Pads</b>                                                                                  |       |                                                                                                                         |            |
| <input type="checkbox"/> 400DI34                                                                                    | —     | Calf Pads, Pair, Standard Integral Foam                                                                                 | NC         |
| <input type="checkbox"/> 400DI35                                                                                    | —     | Calf Pads, Pair, Padded <i>(N/A w/ std. footplate)</i>                                                                  | \$408.00   |
| <b>Angle Adjustable Centermount Calf Pads</b>                                                                       |       |                                                                                                                         |            |
| <input type="checkbox"/> 400DI49                                                                                    | —     | Calf Pads, Pair, Standard Integral Foam                                                                                 | NC         |
| <input type="checkbox"/> 400DI50                                                                                    | —     | Extra Large Single Calf Pad                                                                                             | \$175.00   |
| <b>SEDEO® PRO™ Swing-Away Legrest (Pair) (Footplates, calf supports &amp; options must be selected separately.)</b> |       |                                                                                                                         | <b>\$</b>  |
| <input type="checkbox"/> 400H177                                                                                    | —     | Legrests, Swing-Away, 70°, SEDEO® PRO™ <i>(14" to 21" Leg Length)</i>                                                   | NC         |
| <input type="checkbox"/> 400H178                                                                                    | E0990 | Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ <i>(Includes Padded Calf Pads) (14.5" to 21.5" Leg Length)</i>      | \$408.00   |
| <input type="checkbox"/> 400H179                                                                                    | E1010 | Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>       | \$2,508.00 |
| <input type="checkbox"/> 400HI76                                                                                    | —     | Omit Legrests                                                                                                           | NC         |

## POSITIONING SUPPORTS & ACCESSORIES

| OPTION                                           | HCPC  | DESCRIPTION                                                                                            | RETAIL   |
|--------------------------------------------------|-------|--------------------------------------------------------------------------------------------------------|----------|
| <b>Swing-Away Legrest Footplates</b>             |       |                                                                                                        |          |
| <input type="checkbox"/> 400F369                 | K0040 | Angle Adjustable Flip-Up Footplates                                                                    | \$197.00 |
| <b>Swing-Away Legrest Calf Supports</b>          |       |                                                                                                        |          |
| <input type="checkbox"/> 400DI42                 | —     | Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, depending on chair width, x 9.5" Tall)</i> | \$460.00 |
| <b>Swing-Away Legrest Options</b>                |       |                                                                                                        |          |
| <input type="checkbox"/> 400DI44                 | E0951 | Heel Loops                                                                                             | \$58.00  |
| <input type="checkbox"/> 400KG4                  | —     | Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>                                                | \$234.00 |
| <b>SEDEO® Lower Extremity Swing-Away Support</b> |       |                                                                                                        |          |
| <input type="checkbox"/> 400VD01                 | E1020 | Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right                                             | \$460.00 |
| <input type="checkbox"/> 400VD02                 | E1020 | Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left                                              | \$460.00 |
| <b>Pelvic Positioning Belts</b>                  |       |                                                                                                        |          |
| <input type="checkbox"/> 400A837                 | —     | Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"                                                | NC       |
| <input type="checkbox"/> 400A839                 | —     | Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"                                                | NC       |
| <input type="checkbox"/> 400A840                 | E0978 | Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"                                         | \$88.00  |
| <input type="checkbox"/> 400A845                 | K0108 | Bodypoint – Evoflex Pelvic Stabilizer Belt                                                             | \$203.00 |
| <input type="checkbox"/> 400A846                 | K0108 | Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension                                              | \$224.00 |
| <b>Accessories</b>                               |       |                                                                                                        | \$       |
| <b>Ventilator Option</b>                         |       |                                                                                                        |          |
| <input type="checkbox"/> 400VT12                 | E1029 | Vent Tray Fixed <i>(Pulmonetic LTV Models)</i>                                                         | \$542.00 |
| <b>O2 Tank Option</b>                            |       |                                                                                                        |          |
| <input type="checkbox"/> 400A718                 | E2208 | O2 Tank Holder Mounted to Backrest <i>(D and E Styles)</i>                                             | \$379.00 |
| <b>Wheelchair Tray Table</b>                     |       |                                                                                                        |          |
| <input type="checkbox"/> 4007A09                 | E0950 | Tray, Swing-Away                                                                                       | \$408.00 |
| <b>Medial Thigh Support</b>                      |       |                                                                                                        |          |
| <input type="checkbox"/> 4008A80                 | E0957 | Medial Thigh Support                                                                                   | \$197.00 |
| <input type="checkbox"/> 4008A81                 | E1028 | Pommel Swing-Away Hardware                                                                             | \$292.00 |
| <b>Backpack</b>                                  |       |                                                                                                        |          |
| <input type="checkbox"/> 400A1                   | —     | Black Utility                                                                                          | \$63.00  |
| <b>Misc. Accessories</b>                         |       |                                                                                                        |          |
| <input type="checkbox"/> 4008A59                 | —     | Transfer Bar <i>(Pair)</i>                                                                             | \$408.00 |
| <input type="checkbox"/> 4008A52                 | —     | Ram Mount Phone                                                                                        | \$100.00 |
| <input type="checkbox"/> 4008A53                 | —     | Ram Mount Tablet                                                                                       | \$408.00 |
| <input type="checkbox"/> 4008A54                 | —     | Ram Mount Cup Holder                                                                                   | \$100.00 |
| <input type="checkbox"/> 4008A58                 | —     | Rear View Mirror                                                                                       | \$100.00 |
| <input type="checkbox"/> 400A611                 | —     | Crutch Holder                                                                                          | \$147.00 |
| <input type="checkbox"/> 4008A73                 | —     | Bag for Armrest                                                                                        | \$100.00 |
| <input type="checkbox"/> 4008A75                 | —     | Reflectors <i>(Backrest shroud mounted)</i>                                                            | \$100.00 |
| <input type="checkbox"/> 4009A10                 | —     | Utility Hooks                                                                                          | \$60.00  |

| JOYSTICK / ELECTRONICS           |                |       |                                                                                                                                                                                                                                                                                        |            |
|----------------------------------|----------------|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| OPTION                           | HCPC           | TYPE  | DESCRIPTION                                                                                                                                                                                                                                                                            | RETAIL     |
| Electronics Package              |                |       |                                                                                                                                                                                                                                                                                        | \$         |
| <input type="checkbox"/> 400EL01 | —              | VR-2  | VR-2 or R-net Electronics                                                                                                                                                                                                                                                              | NC         |
| <input type="checkbox"/> 400EL02 | —              | LiNX  | LiNX Electronics <i>(Not available with Recline Only. SitSmarter not available on Tilt Only Configuration)</i>                                                                                                                                                                         | NC         |
| VR-2 Joystick                    |                |       |                                                                                                                                                                                                                                                                                        | \$         |
| <input type="checkbox"/> 400J127 | —              | VR-2  | VR-2 Joystick                                                                                                                                                                                                                                                                          | NC         |
| R-net Joystick (Remote)          |                |       |                                                                                                                                                                                                                                                                                        | \$         |
| <input type="checkbox"/> 400J121 | —              | R-net | Standard LED Joystick <i>(w/o Screen Display)</i>                                                                                                                                                                                                                                      | NC         |
| <input type="checkbox"/> 400J116 | —              | R-net | Color Joystick <i>(2.3" Color LCD Screen Display, CJSM1)</i>                                                                                                                                                                                                                           | NC         |
| <input type="checkbox"/> 400J118 | —              | R-net | High Resolution Color Joystick, (CJSM2+BT)<br><i>3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control</i>                                                                                              | \$1,230.00 |
| <input type="checkbox"/> 400J119 | —              | R-net | Omit Remote Control <i>(Required with arm mounted specialty controls)</i>                                                                                                                                                                                                              | NC         |
| LiNX Joystick                    |                |       |                                                                                                                                                                                                                                                                                        | \$         |
| <input type="checkbox"/> 400J143 | —              | LiNX  | REM110 - Drive Only Remote Joystick <i>(Adaptive Load Compensation)</i>                                                                                                                                                                                                                | NC         |
| <input type="checkbox"/> 400J144 | —              | LiNX  | REM200 Series - LED Joystick <i>(Adaptive Load Compensation, LED icon-based menu, Two Rocker Switches) (REM200 not offered with multi-power. See configurator for full offering.)</i>                                                                                                  | NC         |
| <input type="checkbox"/> 400J145 | —              | LiNX  | REM400 - Touch Screen Color Display Joystick <i>(3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Two Stereo Jacks, Audible Cues, Bluetooth, &amp; Adaptive Load Compensation)</i>                                                                  | \$1,230.00 |
| <input type="checkbox"/> 400J148 | —              | LiNX  | REM400 - Touch Screen Color Display Joystick - N/C for multipower <i>(3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Two Stereo Jacks, Audible Cues, Bluetooth, &amp; Adaptive Load Compensation)</i>                                             | NC         |
| <input type="checkbox"/> 400J146 | —              | LiNX  | REM400 - Touch Screen Color Display Joystick with Toggle Switches <i>(3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Programmable Toggle Switches, Two Stereo Jacks, Audible Cues, Bluetooth, &amp; Adaptive Load Compensation)</i>               | \$1,630.00 |
| <input type="checkbox"/> 400J149 | —              | LiNX  | REM400 - Touch Screen Color Display Joystick with Toggle Switches - Multi-Power <i>(3.5" Color Display, Highly Configurable Touch Interface, Two Multi-Purpose Buttons, Programmable Toggle Switches, Two Stereo Jacks, Audible Cues, Bluetooth, &amp; Adaptive Load Compensation)</i> | \$400.00   |
| <input type="checkbox"/> 400J147 | —              | LiNX  | Omit Remote Control <i>(Required with arm mounted specialty controls)</i>                                                                                                                                                                                                              | NC         |
| Joystick Mounting                |                |       |                                                                                                                                                                                                                                                                                        | \$         |
| <input type="checkbox"/> 400JM34 | —              |       | Fixed on Armrest                                                                                                                                                                                                                                                                       | NC         |
| <input type="checkbox"/> 400JM35 | E1028<br>E1032 |       | Parallel Swing-Away                                                                                                                                                                                                                                                                    | \$345.00   |
| <input type="checkbox"/> 400JM36 | E1028<br>E1032 |       | Link-It Mid Line Mount on Armrest <i>(Swing-Away)</i>                                                                                                                                                                                                                                  | \$991.00   |
| <input type="checkbox"/> 400JM49 | —              |       | Attendant <i>(Available with R-net LED or CJSM1 or LiNX REM100, REM200, or REM400 only. Mounted to Push Handle)</i>                                                                                                                                                                    | NC         |
| Joystick Mount (Location)        |                |       |                                                                                                                                                                                                                                                                                        | \$         |
| <input type="checkbox"/> 400HM16 | —              |       | Right Hand                                                                                                                                                                                                                                                                             | NC         |
| <input type="checkbox"/> 400HM17 | —              |       | Left Hand                                                                                                                                                                                                                                                                              | NC         |
| Display Language Selection       |                |       |                                                                                                                                                                                                                                                                                        | \$         |
| <input type="checkbox"/> 400L1   | —              |       | English                                                                                                                                                                                                                                                                                | NC         |
| <input type="checkbox"/> 400L3   | —              |       | French                                                                                                                                                                                                                                                                                 | NC         |
| <input type="checkbox"/> 400L7   | —              |       | Spanish                                                                                                                                                                                                                                                                                | NC         |

## JOYSTICK / ELECTRONICS

| OPTION | HCPC | TYPE | DESCRIPTION | RETAIL |
|--------|------|------|-------------|--------|
|--------|------|------|-------------|--------|

## Joystick Handle (Optional) \$

|                          |         |       |       |                                              |          |
|--------------------------|---------|-------|-------|----------------------------------------------|----------|
| <input type="checkbox"/> | 400JH13 | E2323 | Both  | Ball Handle                                  | \$147.00 |
| <input type="checkbox"/> | 400JH15 | E2323 | R-net | T-Handle                                     | \$147.00 |
| <input type="checkbox"/> | 400JH24 | E2323 | Both  | Bodypoint – Rubber Dome Joystick Handle      | \$191.00 |
| <input type="checkbox"/> | 400JH25 | E2323 | Both  | Bodypoint – U-Shaped Joystick Handle         | \$191.00 |
| <input type="checkbox"/> | 400JH26 | E2323 | Both  | Bodypoint – Grooved Mushroom Joystick Handle | \$191.00 |
| <input type="checkbox"/> | 400JH30 | E2323 | LiNX  | Chin Cup                                     | \$147.00 |
| <input type="checkbox"/> | 400JH31 | E2323 | LiNX  | Golf Ball                                    | \$147.00 |

## Electronic Programmers \$

|                          |         |   |                              |          |
|--------------------------|---------|---|------------------------------|----------|
| <input type="checkbox"/> | 400PR13 | – | R-net PC Wireless Programmer | \$613.00 |
| <input type="checkbox"/> | 400PR14 | – | LiNX Access Key (LAK)        | \$225.00 |

## Power Modules \$

|                          |         |       |       |                                                                                               |            |
|--------------------------|---------|-------|-------|-----------------------------------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 400ER12 | –     | VR-2  | Power Module, VR-2 90A, Non-Expandable                                                        | NC         |
| <input type="checkbox"/> | 400ER22 | –     | R-net | Power Module, R-net, 90A, Non-Expandable <i>(Only available with Single Power)</i>            | NC         |
| <input type="checkbox"/> | 400ER13 | –     | R-net | Power Module, R-net, 90A, Expandable <i>(Required for Specialty Controls or 2+ Actuators)</i> | \$1,335.00 |
|                          | 400EX11 | E2377 | R-net | Expandable 90A Controller - R-net                                                             | \$706.00   |
|                          | 400EX12 | E2313 | R-net | Harness for 90A Expandable Controller                                                         | \$629.00   |
| <input type="checkbox"/> | 400ER24 | –     | LiNX  | Power Module, LiNX, 120A, Non-Expandable <i>(Only available with Single Power)</i>            | NC         |
| <input type="checkbox"/> | 400ER25 | –     | LiNX  | Power Module, LiNX, 120A, Expandable <i>(Required for Specialty Controls or 2+ Actuators)</i> | \$1,335.00 |
|                          | 400EX19 | E2377 | LiNX  | Expandable 120A Controller - LiNX                                                             | \$706.00   |
|                          | 400EX20 | E2313 | LiNX  | Harness for 120A Expandable Controller - LiNX                                                 | \$629.00   |

## Seating Control \$

## Powered Seating Control

|                          |         |       |                                                  |            |
|--------------------------|---------|-------|--------------------------------------------------|------------|
| <input type="checkbox"/> | 400E040 | –     | Omit Thru-Drive <i>(Requires Switch Control)</i> | NC         |
| <input type="checkbox"/> | 400E037 | E2310 | Thru-drive Control for 1 Actuator                | \$2,240.00 |
| <input type="checkbox"/> | 400E038 | E2311 | Thru-drive Control for 2+ Actuators              | \$4,055.00 |

Switch Control *(Available with Omit Thru-Drive only) (1 Actuator for Tilt or Elevate only)*

|                          |         |   |       |                                                                                   |    |
|--------------------------|---------|---|-------|-----------------------------------------------------------------------------------|----|
| <input type="checkbox"/> | 400E033 | – | R-net | Dual Toggle for 1 Actuator                                                        | NC |
| <input type="checkbox"/> | 400E034 | – | Both  | 2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i> | NC |
| <input type="checkbox"/> | 400E041 | – | R-net | Ctrl+5 Switch Box with Voice Activation for 2 Actuators or more                   | NC |
| <input type="checkbox"/> | 400E042 | – | LiNX  | LiNX 10 Way Switch Box for Actuator Control for 1 Actuator or more                | NC |

Additional Power Seating Control *(Available with Thru-Drive only)*

|                          |         |   |                                                                   |            |
|--------------------------|---------|---|-------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 400QG11 | – | Ctrl+5 Switch Box with Voice Activation                           | \$1,459.00 |
| <input type="checkbox"/> | 400QG12 | – | Ctrl+5 Switch Box with Voice Activation <i>(Multi-power only)</i> | NC         |
| <input type="checkbox"/> | 400QG16 | – | LiNX 10 Way Switch Box for Actuator Control                       | NC         |

## Dual Toggle Mounting Location

|                          |        |   |                     |    |
|--------------------------|--------|---|---------------------|----|
| <input type="checkbox"/> | 400TN1 | – | Toggle Mount, Right | NC |
| <input type="checkbox"/> | 400TN2 | – | Toggle Mount, Left  | NC |

## JOYSTICK / ELECTRONICS

| OPTION | HCPC | TYPE | DESCRIPTION | RETAIL |
|--------|------|------|-------------|--------|
|--------|------|------|-------------|--------|

**Power Seating Control Mounting Location** *(Will swing away if mounted on the same side as the joystick)*

|                          |         |   |                  |          |
|--------------------------|---------|---|------------------|----------|
| <input type="checkbox"/> | 400EY19 | – | Right Armrest    | NC       |
| <input type="checkbox"/> | 400EY20 | – | Left Armrest     | NC       |
| <input type="checkbox"/> | 400EY21 | – | Back Cane, Right | \$175.00 |
| <input type="checkbox"/> | 400EY24 | – | Back Cane, Left  | \$175.00 |

**Ctrl+5 Switch Type**

|                          |         |       |                       |          |
|--------------------------|---------|-------|-----------------------|----------|
| <input type="checkbox"/> | 400EY22 | –     | R-net Push Buttons    | NC       |
| <input type="checkbox"/> | 400EY23 | K0108 | R-net Toggle Switches | \$350.00 |

**Attendant Controls**

|                          |         |       |                                                 |            |
|--------------------------|---------|-------|-------------------------------------------------|------------|
| <input type="checkbox"/> | 4000A05 | E2331 | R-net Dual Control Module for Attendant         | \$1,517.00 |
| <input type="checkbox"/> | 4000A06 | E2331 | R-net Color LCD (CJSM1), Control Attendant Only | \$1,517.00 |
| <input type="checkbox"/> | 4000A11 | E2331 | LiNX Attendant Control Unit (ACU200)            | \$995.00   |

**Attendant Mounting Location**

|                          |        |   |                 |    |
|--------------------------|--------|---|-----------------|----|
| <input type="checkbox"/> | 400DN5 | – | Right Back Cane | NC |
| <input type="checkbox"/> | 400DN6 | – | Left Back Cane  | NC |

**R-net Programming****Assignable Buttons**

*Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.*

|                          |         |   |                            |    |
|--------------------------|---------|---|----------------------------|----|
| <input type="checkbox"/> | 400PM12 | – | Include Assignable Buttons | NC |
|--------------------------|---------|---|----------------------------|----|

**Power Centermount Footplate Extension Option**

|                          |         |       |                                                                                                                                                                                                                |            |
|--------------------------|---------|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 400IE03 | K0108 | Include 7.5" Powered Footplate Extension for Transfer Assist <i>(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation</i> | \$1,459.00 |
|--------------------------|---------|-------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|

**Memory Seating Positions**

*Memory Seating Positions - SitSmarter - allows for positions to be saved and then can be accessed with a single command, moving all related actuators at once to the preset location. Positions can be locked for clinical goals and have others unlocked for easily adjustable positions for comfort and functionality that can be adjusted by the user. Angles can be shown on screen to provide live feedback.*

|                          |         |       |                                                                                                             |    |
|--------------------------|---------|-------|-------------------------------------------------------------------------------------------------------------|----|
| <input type="checkbox"/> | 400PM19 | K0108 | LiNX Memory Seating Positions – SitSmarter <i>(SitSmarter not compatible with Tilt Only configurations)</i> | NC |
|--------------------------|---------|-------|-------------------------------------------------------------------------------------------------------------|----|

**R-net Omni2 Options** *For Expandable Controls, choose all that apply.*

|                          |         |   |                                                                                                                                                         |            |
|--------------------------|---------|---|---------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 4000N14 | – | R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth®<br>Choose Mount Below                                    | \$2,072.00 |
| <input type="checkbox"/> | 4000N15 | – | R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth®<br>Bundled Omni2 Available w/ Switch-It Input Devices | NC         |

**LiNX REM500 Display and IN500 Input Module Options**

|                          |         |   |                                                                                                                                                                                             |            |
|--------------------------|---------|---|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 400AI01 | – | LiNX REM500 Only Input Control Module and 3.5" Touch LCD Display with 2 Stereo Jacks, 2 Buttons, Audible Cues and Bluetooth (Choose Mount Below)                                            | \$2,072.00 |
| <input type="checkbox"/> | 400AI02 | – | LiNX REM500 Bundled Input Control Module and 3.5" Touch LCD Display with 2 Stereo Jacks, 2 Buttons, Audible Cues and Bluetooth <i>(Bundled REM500 available w/ Switch-It Input Devices)</i> | NC         |

**JOYSTICK / ELECTRONICS**

| OPTION                                                                                      | HCPC           | TYPE | DESCRIPTION                                                                              | RETAIL   |
|---------------------------------------------------------------------------------------------|----------------|------|------------------------------------------------------------------------------------------|----------|
| <b>Display Mount</b> <i>If used w/ Standard Joystick, it must be mounted opposite side.</i> |                |      |                                                                                          | \$       |
| <input type="checkbox"/> 400DS14                                                            | E1028<br>E1032 |      | Link-It with Parallel Swing-Away Mount w/Rear Mounted Interface, Right (OMNI2 or REM500) | \$578.00 |
| <input type="checkbox"/> 400DS15                                                            | E1028<br>E1032 |      | Link-It with Parallel Swing-Away Mount w/Rear Mounted Interface, Left (OMNI2 or REM500)  | \$578.00 |

**SPECIALTY INPUT DEVICES**

| OPTION                                                                                                      | HCPC    | DESCRIPTION    | RETAIL                                                                   |            |
|-------------------------------------------------------------------------------------------------------------|---------|----------------|--------------------------------------------------------------------------|------------|
| Hand Control Packages <i>Requires Standard Joystick Omit / Only one hand control device can be ordered.</i> |         |                | \$                                                                       |            |
| Proportional Input Devices                                                                                  |         |                |                                                                          |            |
| <input type="checkbox"/>                                                                                    | 400EU14 | E2312          | Switch-It MicroGuide                                                     | \$3,844.00 |
| <input type="checkbox"/>                                                                                    | 400EU17 | E2328          | Switch-It VersaGuide                                                     | \$2,118.00 |
| <input type="checkbox"/>                                                                                    | 400EU18 | E2328          | Switch-It VersaGuide EZ                                                  | \$2,502.00 |
| Mounting Position <i>(Required)</i>                                                                         |         |                |                                                                          |            |
| <input type="checkbox"/>                                                                                    | 400EM04 | —              | Link-It to the Backrest <i>(Not Available w/ Omit Back Cushion Type)</i> | \$1,109.00 |
| <input type="checkbox"/>                                                                                    | 400EM05 | —              | Link-It to the Fixed Joystick Mount                                      | \$345.00   |
| <input type="checkbox"/>                                                                                    | 400EM06 | E1028<br>E1032 | Link-It to the Parallel Swing-Away Joystick Mount                        | \$578.00   |
| <input type="checkbox"/>                                                                                    | 400EM07 | E1028<br>E1032 | Link-It to the Midline Swing-Away                                        | \$1,225.00 |
| <input type="checkbox"/>                                                                                    | 400EM08 | —              | Omit Mount                                                               | NC         |
| Mounting Platforms <i>(Optional)</i> <i>Not available w/ Link-It to Backrest</i>                            |         |                |                                                                          |            |
| <input type="checkbox"/>                                                                                    | 400MG16 | E0950          | SunDisk Hand Tray <i>(Only available w/ MicroSeries Joysticks)</i>       | \$345.00   |
| <input type="checkbox"/>                                                                                    | 400MG15 | E0950          | Bullet Hand Tray                                                         | \$345.00   |
| Mounting Location <i>(Required)</i>                                                                         |         |                |                                                                          |            |
| <input type="checkbox"/>                                                                                    | 400IM09 | —              | Right                                                                    | NC         |
| <input type="checkbox"/>                                                                                    | 400IM10 | —              | Left                                                                     | NC         |
| Switch Options                                                                                              |         |                |                                                                          |            |
| <input type="checkbox"/>                                                                                    | —       | —              | See Switch Selections <i>(Maximum 2)</i>                                 | —          |

**Head Control Packages** *(Select headrest mount on next page)*

\$

*Only one head control device can be ordered.***Proportional Input Devices**

|                                  |       |                                                          |  |            |
|----------------------------------|-------|----------------------------------------------------------|--|------------|
| <input type="checkbox"/> 400EU25 | E2328 | Switch-It DualPro® Proportional Head Array - Facial Pads |  | \$6,668.00 |
| <input type="checkbox"/> 400EU26 | E2328 | Switch-It DualPro® Proportional Head Array - Spot Pads   |  | \$6,668.00 |
| <input type="checkbox"/> 400EU51 | E2328 | Switch-It Vigo Head Drive                                |  | \$5,648.00 |

**Switched Input Devices**

|                                  |       |                                                    |  |            |
|----------------------------------|-------|----------------------------------------------------|--|------------|
| <input type="checkbox"/> 400EU28 | E2330 | Switch-It 3-Switch Head Array - Short Curved Wings |  | \$5,127.00 |
| <input type="checkbox"/> 400EU30 | E2330 | Switch-It 3-Switch Head Array - Long Curved Wings  |  | \$5,127.00 |

**Switch Options**

|                            |   |                                          |  |   |
|----------------------------|---|------------------------------------------|--|---|
| <input type="checkbox"/> — | — | See Switch Selections <i>(Maximum 2)</i> |  | — |
|----------------------------|---|------------------------------------------|--|---|

## SPECIALTY INPUT DEVICES

| OPTION                                                                | HCPC           | DESCRIPTION                                                                                                                                                                                                                                     | RETAIL     |
|-----------------------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|
| Sip 'N Puff                                                           |                |                                                                                                                                                                                                                                                 | \$         |
| <input type="checkbox"/> 400ID73                                      | E2330          | Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware <i>(Requires headrest mount)</i>                                                                                                                                                   | \$5,827.00 |
| <input type="checkbox"/> 400ID74                                      | K0108          | Switch-It OptiStop with Sip 'N Puff                                                                                                                                                                                                             | \$2,270.00 |
| Headrest Mounts <i>(Applies to Head/Sip 'N Puff control packages)</i> |                |                                                                                                                                                                                                                                                 | \$         |
| <input type="checkbox"/> 400IJ13                                      | E1028<br>E1033 | Whitmyer® Linx Headrest Mount with Quick Release                                                                                                                                                                                                | \$192.00   |
| <input type="checkbox"/> 400IJ14                                      | E1028<br>E1033 | Whitmyer® Linx-2 Headrest Mount with Quick Release                                                                                                                                                                                              | \$234.00   |
| <input type="checkbox"/> 400IJ15                                      | E1028<br>E1033 | Whitmyer® AXYS Headrest Mount with Quick Release                                                                                                                                                                                                | \$123.00   |
| Chin Control Packages                                                 |                |                                                                                                                                                                                                                                                 | \$         |
| <b>Proportional Input Devices</b>                                     |                |                                                                                                                                                                                                                                                 |            |
| <input type="checkbox"/> 400EU34                                      | E2312          | Switch-It MicroGuide <i>(Switches mounted to the headrest only.)</i>                                                                                                                                                                            | \$3,844.00 |
| <input type="checkbox"/> 400EU36                                      | E2328          | Switch-It VersaGuide                                                                                                                                                                                                                            | \$2,118.00 |
| <input type="checkbox"/> 400EU37                                      | E2328          | Switch-It VersaGuide EZ                                                                                                                                                                                                                         | \$2,502.00 |
| <b>Mounting Position</b>                                              |                |                                                                                                                                                                                                                                                 |            |
| <input type="checkbox"/> 400EM09                                      | E1028<br>E1032 | Link-It To Chin Harness, Small (150mm / 6")                                                                                                                                                                                                     | \$1,459.00 |
| <input type="checkbox"/> 400EM10                                      | E1028<br>E1032 | Link-It To Chin Harness, Large (220mm / 8")                                                                                                                                                                                                     | \$1,459.00 |
| <input type="checkbox"/> 400EM11                                      | E1028<br>E1032 | Link-It To Headrest, Swing-Away                                                                                                                                                                                                                 | \$1,459.00 |
| <input type="checkbox"/> 400EM16                                      | —              | Omit Mount                                                                                                                                                                                                                                      | NC         |
| <b>Input Device Mounting Location</b>                                 |                |                                                                                                                                                                                                                                                 |            |
| <input type="checkbox"/> 400IM13                                      | —              | Right                                                                                                                                                                                                                                           | NC         |
| <input type="checkbox"/> 400IM14                                      | —              | Left                                                                                                                                                                                                                                            | NC         |
| <b>Switch Options</b>                                                 |                |                                                                                                                                                                                                                                                 |            |
| <input type="checkbox"/> —                                            | —              | See Switch Selections <i>(Minimum 1 for the Power Chin Boom, Maximum 3)</i>                                                                                                                                                                     | —          |
| Safety Devices                                                        |                |                                                                                                                                                                                                                                                 | \$         |
| <input type="checkbox"/> 400SA06                                      | —              | Switch-It Remote Stop Switch - XLR Connector <i>(Not available with LiNX Electronics)</i>                                                                                                                                                       | \$1,091.00 |
| Switched Driving Packages                                             |                |                                                                                                                                                                                                                                                 | \$         |
| <input type="checkbox"/> 400SW44                                      | E2328          | Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving.<br><i>(Pro Spot proportional switches are required if proportional driving is desired.)<br/>(For non-proportional driving, single switches can be used.)</i> | \$4,048.00 |

## SPECIALTY INPUT DEVICES

## Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.

Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, IN200, IN500, head controls

## Switch ordering instructions:

1. Select switch type and quantity in Section "A" using columns below. (Max 2 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

## Section A – Switch Selections

## Switch / Quantity

| Single Function Switches - Mechanical Switches                                     |                          | Switch 1 |          | \$ | Switch 2                 |          | \$       |
|------------------------------------------------------------------------------------|--------------------------|----------|----------|----|--------------------------|----------|----------|
| 2.5" Single Switch - Black                                                         | <input type="checkbox"/> | 400IR31  | \$227.00 |    | <input type="checkbox"/> | 400IR56  | \$227.00 |
| 2.5" Single Switch - Blue                                                          | <input type="checkbox"/> | 400IR32  | \$227.00 |    | <input type="checkbox"/> | 400IR57  | \$227.00 |
| 2.5" Single Switch - Green                                                         | <input type="checkbox"/> | 400IR33  | \$227.00 |    | <input type="checkbox"/> | 400IR58  | \$227.00 |
| 2.5" Single Switch - Red                                                           | <input type="checkbox"/> | 400IR34  | \$227.00 |    | <input type="checkbox"/> | 400IR59  | \$227.00 |
| 2.5" Single Switch - Yellow                                                        | <input type="checkbox"/> | 400IR35  | \$227.00 |    | <input type="checkbox"/> | 400IR60  | \$227.00 |
| Egg Switch - Black                                                                 | <input type="checkbox"/> | 400IR36  | \$129.00 |    | <input type="checkbox"/> | 400IR61  | \$129.00 |
| Egg Switch - Blue                                                                  | <input type="checkbox"/> | 400IR37  | \$129.00 |    | <input type="checkbox"/> | 400IR62  | \$129.00 |
| Egg Switch - Green                                                                 | <input type="checkbox"/> | 400IR38  | \$129.00 |    | <input type="checkbox"/> | 400IR63  | \$129.00 |
| Egg Switch - Red                                                                   | <input type="checkbox"/> | 400IR39  | \$129.00 |    | <input type="checkbox"/> | 400IR64  | \$129.00 |
| Egg Switch - Yellow                                                                | <input type="checkbox"/> | 400IR40  | \$129.00 |    | <input type="checkbox"/> | 400IR65  | \$129.00 |
| MicroLite (Includes color labels)                                                  | <input type="checkbox"/> | 400IR41  | \$192.00 |    | <input type="checkbox"/> | 400IR66  | \$192.00 |
| Mini Button 1 - Blue (Velcro)                                                      | <input type="checkbox"/> | 400IR42  | \$123.00 |    | <input type="checkbox"/> | 400IR67  | \$123.00 |
| Mini Button 1 - Green (Velcro)                                                     | <input type="checkbox"/> | 400IR43  | \$123.00 |    | <input type="checkbox"/> | 400IR68  | \$123.00 |
| Mini Button 1 - Red (Velcro)                                                       | <input type="checkbox"/> | 400IR44  | \$123.00 |    | <input type="checkbox"/> | 400IR69  | \$123.00 |
| 1" Single Switch - Black                                                           | <input type="checkbox"/> | 400IR45  | \$321.00 |    | <input type="checkbox"/> | 400IR70  | \$321.00 |
| 1" Single Switch - Blue                                                            | <input type="checkbox"/> | 400IR46  | \$321.00 |    | <input type="checkbox"/> | 400IR71  | \$321.00 |
| 1" Single Switch - Green                                                           | <input type="checkbox"/> | 400IR47  | \$321.00 |    | <input type="checkbox"/> | 400IR72  | \$321.00 |
| 1" Single Switch - Orange                                                          | <input type="checkbox"/> | 400IR48  | \$321.00 |    | <input type="checkbox"/> | 400IR73  | \$321.00 |
| 1" Single Switch - Red                                                             | <input type="checkbox"/> | 400IR49  | \$321.00 |    | <input type="checkbox"/> | 400IR74  | \$321.00 |
| 1" Single Switch - Yellow                                                          | <input type="checkbox"/> | 400IR50  | \$321.00 |    | <input type="checkbox"/> | 400IR75  | \$321.00 |
| <b>Electronic Switches</b>                                                         |                          |          |          |    |                          |          |          |
| Flex Switch (Headrest mount only)                                                  | <input type="checkbox"/> | 400IR54  | \$250.00 |    | <input type="checkbox"/> | 400IR79  | \$250.00 |
| <b>ProSpot - Proportional Switches</b>                                             |                          |          |          |    |                          |          |          |
| ProSpot Proportional Switch - Red                                                  | <input type="checkbox"/> | 400IR181 | \$554.00 |    | <input type="checkbox"/> | 400IR186 | \$554.00 |
| ProSpot Proportional Switch - Green                                                | <input type="checkbox"/> | 400IR182 | \$554.00 |    | <input type="checkbox"/> | 400IR187 | \$554.00 |
| ProSpot Proportional Switch - Blue                                                 | <input type="checkbox"/> | 400IR183 | \$554.00 |    | <input type="checkbox"/> | 400IR188 | \$554.00 |
| ProSpot Proportional Switch - Yellow                                               | <input type="checkbox"/> | 400IR184 | \$554.00 |    | <input type="checkbox"/> | 400IR189 | \$554.00 |
| ProSpot Proportional Switch - Orange                                               | <input type="checkbox"/> | 400IR185 | \$554.00 |    | <input type="checkbox"/> | 400IR190 | \$554.00 |
| <b>Section B – Switch Mounting Position</b>                                        |                          |          |          |    |                          |          |          |
| Link-It to the Backrest (Up to 1 per side)                                         | <input type="checkbox"/> | 400EM17  | \$759.00 |    | <input type="checkbox"/> | 400EM23  | \$759.00 |
| Link-It to the Headrest - Gooseneck (Up to 2 per side)                             | <input type="checkbox"/> | 400EM18  | \$759.00 |    | <input type="checkbox"/> | 400EM24  | \$759.00 |
| Link-It to the Fixed Joystick Mount (Up to 2 per side)                             | <input type="checkbox"/> | 400EM19  | \$578.00 |    | <input type="checkbox"/> | 400EM25  | \$578.00 |
| Link-It to the Parallel Swing-Away Joystick Mount (Up to 2 per side)               | <input type="checkbox"/> | 400EM20  | \$578.00 |    | <input type="checkbox"/> | 400EM26  | \$578.00 |
| Link-It to the Seat Rail (Up to 2 per side)<br>(Not avail w/knee & thigh supports) | <input type="checkbox"/> | 400EM21  | \$460.00 |    | <input type="checkbox"/> | 400EM27  | \$460.00 |
| Omit Mount (includes Velcro)                                                       | <input type="checkbox"/> | 400EM22  | NC       |    | <input type="checkbox"/> | 400EM28  | NC       |
| <b>Section C – Switch Mounting Location</b>                                        |                          |          |          |    |                          |          |          |
| Right side of chair                                                                | <input type="checkbox"/> | 400IM17  | NC       |    | <input type="checkbox"/> | 400IM19  | NC       |
| Left side of chair                                                                 | <input type="checkbox"/> | 400IM18  | NC       |    | <input type="checkbox"/> | 400IM20  | NC       |

## Joystick Integrated Switches (Pair)

|                          |         |   |                                                                             |            |
|--------------------------|---------|---|-----------------------------------------------------------------------------|------------|
| <input type="checkbox"/> | 400J003 | – | Proximity Switches on Flex Mount, Pair (Used with VersaGuide/VersaGuide EZ) | \$1,161.00 |
|--------------------------|---------|---|-----------------------------------------------------------------------------|------------|

## JAY® SYNCRA™ SEATING MEASURING GUIDE

### JAY® SYNCRA™ BACK



| Backrest Width       | 16"   | 18"   | 20"   | 22"   |
|----------------------|-------|-------|-------|-------|
| Max User Trunk Width | 13.5" | 15.5" | 17.5" | 19.5" |


#### JAY® Swing-Away Lateral Thoracic Adjustment Range

| Syncra Back Width | Width Range | Height range from bottom of shell to seat pan |
|-------------------|-------------|-----------------------------------------------|
| 16"               | 10" to 24"  | 12" to 19.5"                                  |
| 18"               | 12" to 26"  | 12.5" to 21.5"                                |
| 20"               | 14" to 28"  | 12.5" to 22"                                  |
| 22"               | 16" to 30"  | 13.5" to 23"                                  |

#### Syncra™ Fixed Lateral Support Accessory Range

| Syncra Back Width | Width Range    | Height range from bottom of shell to seat pan | Contour Depth Range |
|-------------------|----------------|-----------------------------------------------|---------------------|
| 16"               | 16.4" to 16.8" | 6" to 8"                                      | 3" to 5"            |
| 18"               | 18.4" to 18.8" | 6" to 8"                                      | 3" to 5"            |
| 20"               | 20.4" to 20.8" | 6" to 8"                                      | 3" to 5"            |
| 22"               | 22.4" to 22.8" | 6" to 8"                                      | 3" to 5"            |

#### Syncra™ Back Shell Height Diagram



### JAY® SYNCRA™ CUSHION



#### JAY® Lateral Thigh Support Adjustment Range

| Syncra Cushion Width | Width Range    |
|----------------------|----------------|
| 16"                  | 8.5" to 16.0"  |
| 17"                  | 9.5" to 17.0"  |
| 18"                  | 10.5" to 18.0" |
| 19"                  | 11.5" to 19.0" |
| 20"                  | 12.5" to 20.0" |
| 22"                  | 14.5" to 22.0" |



Q300 M with SEDEO® PRO™

\$

|                               |                      |
|-------------------------------|----------------------|
| Base Price                    | <input type="text"/> |
| Total of additional equipment | <input type="text"/> |
| Packaging and transport       | <input type="text"/> |
| Total Price                   | <input type="text"/> |

Signature

Send E-Mail