



Quickie QRi
July 2025
EIR5

Mark For:

Submitting for:

☐ Quote

☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

✓ Option # HCPCS Definition Price ✓ Option # HCPCS Definition Price

1. FRAME STYLE & SETUP

Base Model

☐	202WM81	K0005	QRi (265 lbs) Classic Frame Style E1235 Applies with sizes 14" and below	\$2,234.00
--------------------------	---------	-------	--	------------

Transit

☐	202TO1	K0108	Transit	\$374.00
--------------------------	--------	-------	---------	----------

Frame Angle (see Figure 1)

☐	202FA1	75°	NC
☐	202FA8	80°	NC
☐	202FA2	85°	NC

Frame Inset Per Side (See Figure 2)

☐	202FF1	0"	NC
☐	202FF2	1" (N/A w/ 12" wide)	NC

Frame Inset	Footrest Width
0	Subtract 0 from Frame Width
1"	Subtract 2" from Frame Width

Figure 1





Figure 2



Seat Width

Measure from outside of seat tube to outside of opposite seat tube.

☐	202S3	12"	NC
☐	202S4	13"	NC
☐	202S5	14"	NC
☐	202S6	15"	NC
☐	202S7	16"	NC
☐	202S8	17"	NC
☐	202S9	18"	NC
☐	202S10	19"	NC
☐	202S48	E2201 20"	\$583.00

Seat Depth

Measures from front of back post to front edge of seat sling. Seat sling will be 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.


☐	202SD5	12"	NC
☐	202SD6	13"	NC
☐	202SD7	14"	NC
☐	202SD8	15"	NC
☐	202SD9	16"	NC
☐	202SD10	17"	NC
☐	202SD11	18"	NC
☐	202SD12	19"	NC
☐	202SD13	E2203 20"	\$641.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
2. FRONT SEAT HEIGHT SET UP									

Front Seat Height						
See Figure 3						
☐	202SH12	15.0"	NC	☐	145SH20	18.5"
☐	202SH13	15.5"	NC	☐	145SH21	19.0"
☐	202SH14	16.0"	NC	☐	145SH23	19.5"
☐	202SH15	16.5"	NC	☐	145SH24	20.0"
☐	202SH16	17.0"	NC			
☐	202SH17	17.5"	NC			
☐	202SH18	18.0"	NC			

Frame Depth	Max Seat Squeeze
12"/13"	2.5"
14"/15"	3.0"
16"/17"	3.5"
18"/19"	4.0"
20"	4.5"

Figure 3



Caster Wheels				
☐	202CW3	4" Low Profile Polyurethane		NC
☐	202CW44	E2219 4" x 1.25" Semi Pneumatic		\$63.00
☐	202CW64	4" x 1.5" Aluminum Soft Roll		\$216.00
☐	202CW5	5" Low Profile Polyurethane		NC
☐	202CW65	5" x 1" Aluminum Polyurethane		\$129.00
☐	202CW7	E2219 5" x 1.5" Semi Pneumatic		\$100.00
☐	202CW66	5" x 1.5" Aluminum Soft Roll		\$227.00
☐	202CW34	6" Polyurethane		NC
☐	202CW35	E2219 6" x 1.5" Semi Pneumatic		\$100.00
☐	202CW68	6" x 1.5" Aluminum Soft Roll		\$239.00

Caster Fork				
Frog Legs N/A w/ 6" casters				
☐	202CF26		Standard	NC
☐	202CF53	E1015	Frog Legs 0 weight to 100 lbs	\$507.00
☐	202CF54	E1015	Frog Legs 101 lbs weight to 190 lbs	\$507.00
☐	202CF55	E1015	Frog Legs 191 lbs weight to 265 lbs	\$507.00

3. FOOTPLATE OPTIONS

Footrest				
Leg strap not included, see Footplate Option section				
☐	202F308		Aluminum Tubular	NC
☐	202F309		Aluminum Tubular w/ Snap on Cover	NC
☐	202F310	K0040	Aluminum Platform (Ang. Adj.)	\$110.00

Footrest Options				
☐	202DI19	K0038	Leg Strap	\$47.00
☐	202A129	K0038	Extra Leg Strap	\$47.00
☐	202DI06	K0108	Positioning Plates for Platform Footrest	\$58.00
☐	202A500		Gel Cover Hanger/Frame	\$123.00

Adaptive Platform Spacer				
N/A with Tubular Footrest				
☐	202AD1		Adaptive Platform 2" Spacer	\$110.00
☐	202AD2		Adaptive Platform 4" Spacer	\$110.00
☐	202AD3		Adaptive Platform 6" Spacer	\$110.00

Tubular and Platform Footrest Caster/Fork Availability

Lower Leg Length Range

Frame Angle	Footrest Type	Range	FSH Range 15"-17"		FSH Range 17.5"-20"	
			Inset	0"	1"	0"
75°	Platform	Min		12.3	12.3	14.3
		Max		16.0	16.0	18.0
	Tubular	Min		12.3	12.3	14.3
		Max		15.5	15.5	17.5
	Long Platform and Tubular	Min		12.3	13.9	14.3
		Max		17.5	17.5	19.5
80°	Platform	Min		12.2	12.1	14.2
		Max		16.1	16.1	18.1
	Tubular	Min		12.2	12.2	14.2
		Max		15.6	15.6	17.6
	Long Platform and Tubular	Min		12.2	14.2	14.2
		Max		17.8	17.8	19.8
85°	Platform	Min		12.2	12.5	14.2
		Max		16.1	16.1	18.1
	Tubular	Min		12.2	12.5	14.2
		Max		15.5	15.5	17.5
	Long Platform and Tubular	Min		12.2	14.5	14.2
		Max		17.7	17.7	19.7

			Inset	
			0"	1"
75° Frame Angle	4	Standard Fork		
		Frog Legs		
	5	Standard Fork		
		Frog Legs		
	6	Standard Fork		
		Frog Legs	N/A	N/A
80° Frame Angle	4	Standard Fork		
		Frog Legs		
	5	Standard Fork		
		Frog Legs	N/A	
	6	Standard Fork	N/A w/1.5	
		Frog Legs	N/A	N/A
85° Frame Angle	4	Standard Fork		
		Frog Legs	N/A	N/A w/1.5
	5	Standard Fork	N/A w/1.5	
		Frog Legs	N/A	N/A
	6	Standard Fork	N/A	N/A
		Frog Legs	N/A	N/A

N/A w/ 1.5 = N/A with the 1.5" wide caster. Blank = available.
No Go's may still apply that are not listed based on actual configuration of the chair.

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
4. REAR SEAT HEIGHT SET UP									

Rear Seat Height				
<i>14" and 14.5" only available w/ 22" rear wheels</i>				
☐	202RH13	14.0"		NC
☐	202RH14	14.5"		NC
☐	202RH15	15.0"		NC
☐	202RH16	15.5"		NC
☐	202RH17	16.0"		NC
☐	202RH18	16.5"		NC
☐	202RH19	17.0"		NC
☐	202RH26	17.5"		NC
☐	202RH20	18.0"		NC
☐	202RH27	18.5"		NC
☐	202RH21	19.0"		NC
☐	202RH37	19.5"		NC
☐	202RH22	20.0"		NC

Rear Wheel Size				
<i>Ultra Lightweight Wheels only available with aluminum anodized handrim</i>				
<i>Contact Built 4 Me for Spinergy color options</i>				
☐	202RW8	22" Mag Wheel		NC
☐	202RW1	24" Mag Wheel		NC
☐	202RW91	22" Lite Spoke Wheel	\$234.00	
☐	202RW92	24" Lite Spoke Wheel		NC
☐	202RW93	25" Lite Spoke Wheel		NC
☐	202RW94	26" Lite Spoke Wheel	\$234.00	
☐	202RW06	22" Ultra Lightweight Spoke Wheel 2.0	\$292.00	
☐	202RW10	24" Ultra Lightweight Spoke Wheel 2.0	\$292.00	
☐	202RW16	25" Ultra Lightweight Spoke Wheel 2.0	\$292.00	
☐	202RW28	22" Spinergy SPOX Black	\$985.00	
☐	202RW95	24" Spinergy SPOX Black	\$985.00	
☐	202RW96	25" Spinergy SPOX Black	\$985.00	
☐	202RW97	26" Spinergy SPOX Black	\$985.00	
☐	202RW32	22" Spinergy LX Black	\$1,237.00	
☐	202RW98	24" Spinergy LX Black	\$1,237.00	
☐	202RW99	25" Spinergy LX Black	\$1,237.00	
☐	202RW40	26" Spinergy LX Black	\$1,237.00	

Rear Wheel Tire				
<i>See Availability Matrix</i>				
☐	202RT1	Pneumatic		NC
☐	202RT2	E2211 E2213	Pneumatic w/ Airless Insert	\$100.00
☐	202RT3		Full Poly	NC
☐	202RT6	E2211 E2212	High-Pressure Clincher	\$110.00
☐	202RT27	E2211 E2212	Knobby V- Trak	\$110.00
☐	202RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	\$164.00
☐	202RT38	E2218	SHOX Performance Poly	\$110.00

Short Axle Plate/Wheel Size Matrix				
Rear Seat to Floor Height	22"	24"	25"	26"
	14.0			
	14.5	15.0	15.0	16.0
	16.0	15.5	15.5	17.5
	16.5	17.0	17.0	18.0
	17.0	17.5	17.5	18.5
	17.5	18.0	18.0	19.0

Rear Tire & Wheel Matrix				
Tire	Wheel	Mag	Lite Spoke	Ultra Lightweight
Pneumatic		22,24	22,24,26	
Pneumatic w/ Airless Insert		22,24	22,24	
Full Poly		22,24	22,24	
High-Pressure Clincher			22,24,25,26	22,24,25
Knobby V Trak by Kenda				24
Schwalbe Marathon			24,25,26	24,25,26
SHOX Performance Poly	24	24,25		24,25
22" = 489 24" = 540 25" = 559 26" = 590				

Handrim				
<i>Plastic not available on all wheels in all sizes. Natural Fits are N/A w/ 26"</i>				
☐	202HR1	Aluminum Anodized		NC
☐	202HR4	Plastic		\$129.00
☐	202HR22	Natural Fit - Standard Grip		\$549.00
☐	202HR23	Natural Fit - Super Grip		\$611.00

Center of Gravity Preset				
<i>Measure from front of back post to center of rear axle</i>				
☐	202CC2	1.0"		NC
☐	202CC3	1.5"		NC
☐	202CC4	2.0"		NC
☐	202CC5	2.5"		NC
☐	202CC6	3.0"		NC

Camber				
☐	202CI1	0°		NC
☐	202CI2	3°		NC

Axle				
☐	202AX1	Quick Release Steel		NC
☐	202AX13	Quad Quick Release		\$147.00
<i>Only available with Mag wheels</i>				

Axle Plates				
<i>See Short Axle Plate Matrix for RSH Availability</i>				
☐	202AL40	Standard Short (COG 1" to 3")		NC
☐	202AL43	Standard Long (COG 1" to 3")		NC

Wheel Lock				
<i>Not all wheel locks available in all seat depths, wheel sizes and COG settings</i>				
☐	202WL1	High Mount Push		NC
☐	202WL2	High Mount Pull		NC
☐	202WL54	Ergo Scissor		NC
☐	202WL30	Compact		\$164.00
☐	202WL10	Omit Wheel Lock		NC
☐	202WL12	Do Not Mount		NC

Extension Handle				
<i>Only available w/ High Mount Push or Pull Wheel Locks</i>				
☐	202EH7	E0961	Extension handles (Pair)	\$88.00

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

5. ARMREST

Armrest - Padded Swing Away

Preset to 10". Some limitations may exist with integral and folding push handles Range 9"-13".

☐	202AA18	E0973	Padded Swing-Away - Standard	\$263.00
--------------------------	---------	-------	------------------------------	----------

Armrest - Height Adjustable

☐	202AA11	E0973	Single Post Ht. Adjust. Std (8.5"-13.5")	\$333.00
☐	202AA29	E0973	Single Post Ht. Adjust. Low (7"-11")	\$333.00

Armpad - For Height Adjustable

☐	202AM23		Desk Length - Classic (10")	NC
☐	202AM24		Full Length - Classic (14")	NC

Omit Armrest

☐	202AA14		Omit Armrest	NC
--------------------------	---------	--	--------------	----

Armrest Accessories

☐	2028A05		Gel Cover Swing Away Arm	\$175.00
--------------------------	---------	--	--------------------------	----------

6. BACKREST

Back Type

☐	202BT82		Folding Lock Down Angle Adj	NC
☐	202BT3		Non-Folding Backrest	NC

Push Handles

☐	202PH6		Integral Push Handles	NC
☐	202PH5		Fold-Down Push Handles N/A w/ 10-14" Back Height or Standard Upholstery.	\$116.00
☐	202PH2		Omit Push Handles	NC

Height Adjustable Backrest

☐	202B110		10" - 14" Back Height	NC
☐	202B127		12" - 16" Back Height	NC
☐	202B111		14" - 18" Back Height	NC
☐	202B137		16" - 20" Back Height	NC

Back Upholstery

☐	202BU1		Standard Upholstery	NC
☐	202BU70	E2611	Backrest Tension Adj Ballistic Nylon *reference 186BU70 for billing	\$367.00
☐	202BU71	E2611	Backrest Tension Adj 3DX Vented *reference 186BU71 for billing	\$387.00
☐	202BU14		Omit Upholstery (Use w/ JAY Backs)	-\$15.00

JAY J3 Back

☐	202JB01		JAY J3 Back	See OF
--------------------------	---------	--	-------------	--------

Military Service Patches

Only available on tension adjustable upholstery and Jay 3 Back

☐	202VU1		US Marines Patch	\$116.00
☐	202VU2		US Army Patch	\$116.00
☐	202VU3		US Navy Patch	\$116.00
☐	202VU4		US Air Force Patch	\$116.00
☐	202VU12		US Coast Guard	\$116.00
☐	202VU5		Purple Heart	\$116.00
☐	202VU8		Purple Heart with Marines Patch	\$181.00
☐	202VU9		Purple Heart with Army Patch	\$181.00
☐	202VU10		Purple Heart with Navy Patch	\$181.00
☐	202VU11		Purple Heart with Air Force Patch	\$181.00
☐	202VU13		Purple Heart with Coast Guard	\$181.00

7. SEATING

Seat Sling

☐	202SL1		Standard Seat Sling	NC
☐	202SL6		Omit Seat Sling	-\$15.00

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	202JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width:	Depth:	

✓	Option #	HCP	CS	Definition	Price	✓	Option #	HCP	CS	Definition	Price
---	----------	-----	----	------------	-------	---	----------	-----	----	------------	-------

8. Key Options & Accessories

Anti Tips				
☐	202AT1	E0971	Anti Tips (Pair)	\$175.00
☐	202AT3		Omit Rear Anti Tip	NC
Side Guards				
☐	202SG33	K0108	Plastic - Standard (8.5" High)	\$169.00
☐	202SG29	K0108	Plastic - Low (6.5" High)	\$169.00
☐	202SG35	K0108	Aluminum <i>N/A w/ height adj arms, push/pull wheel locks only available in black</i>	\$286.00
☐	202SG36	K0108	Aluminum w/ Fender <i>N/A w/ height adj arms, push/pull wheel locks only available in black</i>	\$321.00

Pack Accessories				
☐	202A1		Quickie Backpack Black	\$88.00
☐	202A2		Quickie Kids Backpack Black	\$70.00
☐	202A8		Seat Pouch Black	\$58.00
Positioning Belt				
☐	202A141	E0978	1.5" Auto Buckle Universal w/clamps	\$82.00
☐	202A120		Positioning Belt Slide Clamps	\$58.00
Spoke Guards				
☐	202A302	K0065	Quickie Clear Spoke Guard	\$134.00
☐	202A23	K0065	Quickie Black Spoke Guard	\$134.00
Miscellaneous				
☐	202DI07		Neoprene Impact Guards	\$105.00
☐	202A2100	E2207	Crutch Holder <i>Only available with Integral or Fold Down push handles</i>	\$227.00
☐	202A85		Bolt-On Push Handles	\$134.00
☐	202A83		Tool Kit	\$63.00

9. COLOR

NEW!	Standard Frame Colors			
☐	202C214		Black (Gloss)	NC
☐	202C215		Black (Matte)	NC
☐	202C211		Celestial Sapphire	NC
☐	202C212		Desert Sage	NC
☐	202C213		Desert Tan	NC
☐	202C216		Metallic Shadow	NC
☐	202C217		Midnight Blue	NC
☐	202C218		Mystic Purple	NC
☐	202C219		Petal Pink	NC
☐	202C220		Pink Power	NC
☐	202C221		Red Pepper	NC
☐	202C222		Royal Blue	NC
☐	202C223		Salute Green	NC
☐	202C224		Scarlet Noir	NC
☐	202C225		Silver Ore	NC
☐	202C226		Solar Orange	NC
☐	202C227		Starlight White	NC
☐	202C228		Storm Gray	NC
☐	202C229		Turbo Green	NC
☐	202C230		Velocity Yellow	NC

Personalized Frame Color				
☐	202C112		Desert Camouflage	\$460.00
☐	202C43		Mossy Oak Camouflage	\$460.00
☐	202C42		Stars & Stripes	\$460.00
☐	202C108		Zebra (<i>Choose above base color</i>) Std with White unless above base paint is chosen N/A w/ Black, Matte Black)	\$460.00

Additional Options Through Built for Me

To inquire about any unique personalized Quickie products, contact: Phone: (800) 333-4000, Fax: (800) 333-9011



Sunrise Medical (US) LLC
12002 Volunteer Blvd. Mount Juliet, TN 37122