



ORDER FORM

Horizon Stander

EILK02

Robust and durable 3-in-1 standing frame for both young children and adult users.

Customer			
Account			
Account #		Date	
☐ Quote	☐ Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Customer Service: 800-333-4000 | Fax: 800-300-7502 | Email: customseating@sunmed.com | www.sunrisemedical.com

Horizon Stander

Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 220 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Frame			\$

Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals. (Knee supports and Footplates separate)

☐ 141FO5	Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E0641	\$1,948
☐ 141FO6	Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E0641	\$2,030
☐ 141FO7	Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Vinyl	E0641	\$1,948
☐ 141FO8	Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Vinyl	E0641	\$2,030

Base			\$
------	--	--	----

Selection Required. MUST match Frame size. Includes worldwide adapters.

☐ 141NB110	Size 2 Horizon Base/Chassis		\$3,599
☐ 141NB111	Size 3 Horizon Base/Chassis		\$4,491

Headrest			\$
----------	--	--	----

Head Support Assembly selection required if Head Lateral Support/Extension Bracket is selected.

☐ 141HE3	Size 2 and 3 Head Support Assembly - Black		\$293
☐ 141HE4	Size 2 and 3 Head Lateral Supports - Black		\$253
☐ 141HE63	Size 2 and 3 Head Support Assembly - Vinyl		\$293
☐ 141HE64	Size 2 and 3 Head Lateral Supports - Vinyl		\$253
☐ 141JG12	Headrest 4" Horizontal Extension Bracket		\$141
☐ 141Y029	Size 1-3 Contoured Head Support		\$641

Hip Laterals			\$
--------------	--	--	----

Selection required if Hip Belt is not selected.

☐ 141KC7	Size 2 Horizon Hip Laterals - Black (Pair)		\$375
☐ 141KC8	Size 3 Horizon Hip Laterals - Black (Pair)		\$391
☐ 141C29	Size 2 Horizon Hip Laterals - Vinyl (Pair)		\$375
☐ 141C30	Size 3 Horizon Hip Laterals - Vinyl (Pair)		\$391

Hip Belts			\$
-----------	--	--	----

Selection required if Hip Laterals are not selected. Multiple Selections can be made.

☐ 141A856	Size 2 Supine Belt - Black		\$296
☐ 141A857	Size 3 Supine Belt - Black		\$296
☐ 141A883	Size 2 Supine Belt - Vinyl		\$296
☐ 141A884	Size 3 Supine Belt - Vinyl		\$296

Knee Supports			\$
---------------	--	--	----

Selection of Standard Knee Assembly Required.

Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Knee Bracket has 5" - 7" adjustment.

Long Knee Bracket has 7" - 9" adjustment.

☐ 141XX16	Size 2 and 3 Standard Knee Assembly - Black (Pair)		\$526
☐ 141XX17	Size 2 and 3 Standard Knee Assembly - Vinyl (Pair)		\$526
☐ 141JG14	Medium Knee Adjustment Bracket (Pair)		\$121
☐ 141JG15	Long Knee Adjustment Bracket (Pair)		\$124

Lower Leg Positioning			\$
-----------------------	--	--	----

Selection of Footplate Required.

☐ 1418A107	Small Sandals with Straps (Pair)		\$278
☐ 1418A108	Medium Sandals with Straps (Pair)		\$304
☐ 1418A109	Large Sandals with Straps (Pair)		\$327
☐ 141F53	Size 2 and 3 Flip Up Footplate Assembly (Pair)		\$831
☐ 141F65	Size 1-3 3D Footplate Assembly (Pair) - One size fits all		\$1,359

Upper Extremity Support Tray			\$
------------------------------	--	--	----

Must match Horizon Stander Size

☐ 141A209	Size 2 and 3 Supine Tray		\$441
☐ 141A210	Size 2 and 3 Prone Tray		\$583

Technical Specifications - Imperial

Frame Sizing Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	50.4" x 25.2" x 31.7"	50.4" x 25.2" x 31.7"
Weight of Product	110.6 lbs	130 lbs
Maximum User Weight	176 lbs	220 lbs
User Standing Height (Min / Max)	55.1" x 70.9"	61" / 70.9"
Shoulder Height - Footplate to Top of Pad (Min / Max)	39.4" / 57"	39.4" / 57"
Hip Pad Height - Footplate to Top of Pad (Min / Max)	29.5" / 37.4"	29.5" / 45.3"
Hip Pad Width (Min / Max)	8.3" / 15.3"	12.2" / 18.5"
Knee Height - Footplate to Middle of Pad (Min / Max)	15.7" / 21.6"	15.7" / 25.6"
Width Between Chest Laterals (Min / Max)	9.4" / 14.9"	13.4" / 18.1"
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	20.5" x 23.2"	20.5" x 23.2"
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	7"	7"

Technical Specifications - Metric

Frame Sizing Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	1280 mm x 640 mm x 805 mm	1280 mm x 640 mm x 805 mm
Weight of Product	50 kgs	95 kgs
Maximum User Weight	80 kgs	100 kgs
User Standing Height (Min / Max)	1400 mm / 1800 mm	1550 mm / 1800 mm
Shoulder Height - Footplate to Top of Pad (Min / Max)	1000 mm / 1450 mm	1000 mm / 1450 mm
Hip Pad Height - Footplate to Top of Pad (Min / Max)	750 mm / 950 mm	750 mm / 1150 mm
Hip Pad Width (Min / Max)	210 mm / 390 mm	310 mm / 470 mm
Knee Height - Footplate to Middle of Pad (Min / Max)	400 mm / 550 mm	400 mm / 650 mm
Width Between Chest Laterals (Min / Max)	240 mm / 380 mm	340 mm / 460 mm
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	520 mm x 590 mm	520 mm x 590 mm
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	180 mm	180 mm

Sunrise Medical (US) LLC
12002 Volunteer Blvd
Mount Juliet, TN 37122



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129848 Rev. I © 10/2025

