



JAY® GO CUSHION

E2603 / E2604

January 2025

Mark For:

Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for:

Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

JAY® GO Cushion - E2603 / E2604

JAY GO Cushion

✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JGOXXYY	E2603	JAY GO Cushion 14" - 21" wide	\$240
<i>Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)</i>				
<input type="checkbox"/>	JGOXXYY	E2604	JAY GO Cushion 22" - 24" wide	\$361
<i>Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)</i>				

JAY GO Size Matrix

		WIDTH													
		12	13	14	15	16	17	18	19	20	21	22	23	24	
DEPTH	12	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	13	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
	14	N/A	N/A						N/A	N/A	N/A	N/A	N/A	N/A	
	15	N/A	N/A								N/A	N/A	N/A	N/A	
	16	N/A	N/A								N/A	N/A	N/A	N/A	
	17	N/A	N/A										N/A	N/A	
	18	N/A	N/A											N/A	N/A
	19	N/A	N/A												
	20	N/A	N/A												
	21	N/A	N/A	N/A											
	22	N/A	N/A	N/A	N/A										
	23	N/A	N/A	N/A	N/A	N/A									
24	N/A	N/A	N/A	N/A	N/A	N/A									

	=	Standard Offering
	=	Offered Via Jay Your Way
	=	Not available

JAY GO Measuring Guide

Measurements made from completed cushion with covers.

Width

X = Width

Note: Measurement made from outside edge to outside edge of cushion.

Ex JGO1618

X = 16"


Depth

Y = Depth

Note: Measurement made from front edge of rear cutout and front lower edge of cushion.

Ex JGO1618

Y = 18"



JAY GO Cushion

✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JGO1414	E2603	JAY GO Cushion 14" x 14"	\$240
<input type="checkbox"/>	JGO1415	E2603	JAY GO Cushion 14" x 15"	\$240
<input type="checkbox"/>	JGO1416	E2603	JAY GO Cushion 14" x 16"	\$240
<input type="checkbox"/>	JGO1417	E2603	JAY GO Cushion 14" x 17"	\$240
<input type="checkbox"/>	JGO1418	E2603	JAY GO Cushion 14" x 18"	\$240
<input type="checkbox"/>	JGO1419	E2603	JAY GO Cushion 14" x 19"	\$240
<input type="checkbox"/>	JGO1420	E2603	JAY GO Cushion 14" x 20"	\$240
<input type="checkbox"/>	JGO1514	E2603	JAY GO Cushion 15" x 14"	\$240
<input type="checkbox"/>	JGO1515	E2603	JAY GO Cushion 15" x 15"	\$240
<input type="checkbox"/>	JGO1516	E2603	JAY GO Cushion 15" x 16"	\$240
<input type="checkbox"/>	JGO1517	E2603	JAY GO Cushion 15" x 17"	\$240
<input type="checkbox"/>	JGO1518	E2603	JAY GO Cushion 15" x 18"	\$240
<input type="checkbox"/>	JGO1519	E2603	JAY GO Cushion 15" x 19"	\$240
<input type="checkbox"/>	JGO1520	E2603	JAY GO Cushion 15" x 20"	\$240
<input type="checkbox"/>	JGO1521	E2603	JAY GO Cushion 15" x 21"	\$240
<input type="checkbox"/>	JGO1614	E2603	JAY GO Cushion 16" x 14"	\$240
<input type="checkbox"/>	JGO1615	E2603	JAY GO Cushion 16" x 15"	\$240
<input type="checkbox"/>	JGO1616	E2603	JAY GO Cushion 16" x 16"	\$240
<input type="checkbox"/>	JGO1617	E2603	JAY GO Cushion 16" x 17"	\$240
<input type="checkbox"/>	JGO1618	E2603	JAY GO Cushion 16" x 18"	\$240
<input type="checkbox"/>	JGO1619	E2603	JAY GO Cushion 16" x 19"	\$240
<input type="checkbox"/>	JGO1620	E2603	JAY GO Cushion 16" x 20"	\$240
<input type="checkbox"/>	JGO1621	E2603	JAY GO Cushion 16" x 21"	\$240
<input type="checkbox"/>	JGO1622	E2603	JAY GO Cushion 16" x 22"	\$240
<input type="checkbox"/>	JGO1714	E2603	JAY GO Cushion 17" x 14"	\$240
<input type="checkbox"/>	JGO1715	E2603	JAY GO Cushion 17" x 15"	\$240
<input type="checkbox"/>	JGO1716	E2603	JAY GO Cushion 17" x 16"	\$240
<input type="checkbox"/>	JGO1717	E2603	JAY GO Cushion 17" x 17"	\$240
<input type="checkbox"/>	JGO1718	E2603	JAY GO Cushion 17" x 18"	\$240
<input type="checkbox"/>	JGO1719	E2603	JAY GO Cushion 17" x 19"	\$240
<input type="checkbox"/>	JGO1720	E2603	JAY GO Cushion 17" x 20"	\$240
<input type="checkbox"/>	JGO1721	E2603	JAY GO Cushion 17" x 21"	\$240
<input type="checkbox"/>	JGO1722	E2603	JAY GO Cushion 17" x 22"	\$240
<input type="checkbox"/>	JGO1723	E2603	JAY GO Cushion 17" x 23"	\$240
<input type="checkbox"/>	JGO1814	E2603	JAY GO Cushion 18" x 14"	\$240
<input type="checkbox"/>	JGO1815	E2603	JAY GO Cushion 18" x 15"	\$240
<input type="checkbox"/>	JGO1816	E2603	JAY GO Cushion 18" x 16"	\$240
<input type="checkbox"/>	JGO1817	E2603	JAY GO Cushion 18" x 17"	\$240
<input type="checkbox"/>	JGO1818	E2603	JAY GO Cushion 18" x 18"	\$240
<input type="checkbox"/>	JGO1819	E2603	JAY GO Cushion 18" x 19"	\$240
<input type="checkbox"/>	JGO1820	E2603	JAY GO Cushion 18" x 20"	\$240
<input type="checkbox"/>	JGO1821	E2603	JAY GO Cushion 18" x 21"	\$240
<input type="checkbox"/>	JGO1822	E2603	JAY GO Cushion 18" x 22"	\$240
<input type="checkbox"/>	JGO1823	E2603	JAY GO Cushion 18" x 23"	\$240
<input type="checkbox"/>	JGO1824	E2603	JAY GO Cushion 18" x 24"	\$240
<input type="checkbox"/>	JGO1915	E2603	JAY GO Cushion 19" x 15"	\$240
<input type="checkbox"/>	JGO1916	E2603	JAY GO Cushion 19" x 16"	\$240
<input type="checkbox"/>	JGO1917	E2603	JAY GO Cushion 19" x 17"	\$240

✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JGO1918	E2603	JAY GO Cushion 19" x 18"	\$240
<input type="checkbox"/>	JGO1919	E2603	JAY GO Cushion 19" x 19"	\$240
<input type="checkbox"/>	JGO1920	E2603	JAY GO Cushion 19" x 20"	\$240
<input type="checkbox"/>	JGO1921	E2603	JAY GO Cushion 19" x 21"	\$240
<input type="checkbox"/>	JGO1922	E2603	JAY GO Cushion 19" x 22"	\$240
<input type="checkbox"/>	JGO1923	E2603	JAY GO Cushion 19" x 23"	\$240
<input type="checkbox"/>	JGO1924	E2603	JAY GO Cushion 19" x 24"	\$240
<input type="checkbox"/>	JGO2015	E2603	JAY GO Cushion 20" x 15"	\$240
<input type="checkbox"/>	JGO2016	E2603	JAY GO Cushion 20" x 16"	\$240
<input type="checkbox"/>	JGO2017	E2603	JAY GO Cushion 20" x 17"	\$240
<input type="checkbox"/>	JGO2018	E2603	JAY GO Cushion 20" x 18"	\$240
<input type="checkbox"/>	JGO2019	E2603	JAY GO Cushion 20" x 19"	\$240
<input type="checkbox"/>	JGO2020	E2603	JAY GO Cushion 20" x 20"	\$240
<input type="checkbox"/>	JGO2021	E2603	JAY GO Cushion 20" x 21"	\$240
<input type="checkbox"/>	JGO2022	E2603	JAY GO Cushion 20" x 22"	\$240
<input type="checkbox"/>	JGO2023	E2603	JAY GO Cushion 20" x 23"	\$240
<input type="checkbox"/>	JGO2024	E2603	JAY GO Cushion 20" x 24"	\$240
<input type="checkbox"/>	JGO2117	E2603	JAY GO Cushion 21" x 17"	\$240
<input type="checkbox"/>	JGO2118	E2603	JAY GO Cushion 21" x 18"	\$240
<input type="checkbox"/>	JGO2119	E2603	JAY GO Cushion 21" x 19"	\$240
<input type="checkbox"/>	JGO2120	E2603	JAY GO Cushion 21" x 20"	\$240
<input type="checkbox"/>	JGO2121	E2603	JAY GO Cushion 21" x 21"	\$240
<input type="checkbox"/>	JGO2122	E2603	JAY GO Cushion 21" x 22"	\$240
<input type="checkbox"/>	JGO2123	E2603	JAY GO Cushion 21" x 23"	\$240
<input type="checkbox"/>	JGO2124	E2603	JAY GO Cushion 21" x 24"	\$240

22" & Greater in Width

<input type="checkbox"/>	JGO2217	E2604	JAY GO Cushion 22" x 17"	\$361
<input type="checkbox"/>	JGO2218	E2604	JAY GO Cushion 22" x 18"	\$361
<input type="checkbox"/>	JGO2219	E2604	JAY GO Cushion 22" x 19"	\$361
<input type="checkbox"/>	JGO2220	E2604	JAY GO Cushion 22" x 20"	\$361
<input type="checkbox"/>	JGO2221	E2604	JAY GO Cushion 22" x 21"	\$361
<input type="checkbox"/>	JGO2222	E2604	JAY GO Cushion 22" x 22"	\$361
<input type="checkbox"/>	JGO2223	E2604	JAY GO Cushion 22" x 23"	\$361
<input type="checkbox"/>	JGO2224	E2604	JAY GO Cushion 22" x 24"	\$361
<input type="checkbox"/>	JGO2319	E2604	JAY GO Cushion 23" x 19"	\$361
<input type="checkbox"/>	JGO2320	E2604	JAY GO Cushion 23" x 20"	\$361
<input type="checkbox"/>	JGO2321	E2604	JAY GO Cushion 23" x 21"	\$361
<input type="checkbox"/>	JGO2322	E2604	JAY GO Cushion 23" x 22"	\$361
<input type="checkbox"/>	JGO2323	E2604	JAY GO Cushion 23" x 23"	\$361
<input type="checkbox"/>	JGO2324	E2604	JAY GO Cushion 23" x 24"	\$361
<input type="checkbox"/>	JGO2419	E2604	JAY GO Cushion 24" x 19"	\$361
<input type="checkbox"/>	JGO2420	E2604	JAY GO Cushion 24" x 20"	\$361
<input type="checkbox"/>	JGO2421	E2604	JAY GO Cushion 24" x 21"	\$361
<input type="checkbox"/>	JGO2422	E2604	JAY GO Cushion 24" x 22"	\$361
<input type="checkbox"/>	JGO2423	E2604	JAY GO Cushion 24" x 23"	\$361
<input type="checkbox"/>	JGO2424	E2604	JAY GO Cushion 24" x 24"	\$361

JAY Your Way - Modifications


*Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.
Please refer to the "JAY YOUR WAY" product catalog for complete listing of modifications*

Notch Adjustments

- | | | | |
|--------------------------|----|----------------------|----|
| <input type="checkbox"/> | NR | Notch for Seat Rails | NC |
|--------------------------|----|----------------------|----|
- 1.25" from bottom and 1.25" from side of cushion removed.



- | | | | |
|--------------------------|----|------------------------------|----|
| <input type="checkbox"/> | NP | Notch for Posts (Back Canes) | NC |
|--------------------------|----|------------------------------|----|
- From Back of Cushion: 1" 2" 3"
- 1" from side of cushion removed. No Modification to inner or outer cover.



Anterior or Posterior Slope

- | | | | |
|--------------------------|-----|-----------------|------|
| <input type="checkbox"/> | WGP | Anterior Slope | \$29 |
| <input type="checkbox"/> | WGA | Posterior Slope | \$29 |
- Height Increase:
- Foam wedge glued underneath cushion, cover sewn to match.

Adjustments

- | | | | |
|--------------------------|------|------------------------------|----|
| <input type="checkbox"/> | LLCR | Leg Length Cut - Right | NC |
| <input type="checkbox"/> | LLCL | Leg Length Cut - Left | NC |
| <input type="checkbox"/> | CCT | Contracture Cut (1" x 45°) | NC |
| <input type="checkbox"/> | LTC | Leg Trough Cut (1" Standard) | NC |
| <input type="checkbox"/> | MSC | Medial Support Cut | NC |

Cover Options

- | | | | |
|--------------------------|----|-------------|------|
| <input type="checkbox"/> | CC | Cover Color | \$29 |
|--------------------------|----|-------------|------|
- Color:

Embroidery

- | | | | |
|--------------------------|------------|-------------------------------|------|
| <input type="checkbox"/> | EMBROIDERY | Embroidery (32 Character Max) | \$26 |
|--------------------------|------------|-------------------------------|------|
- Font: Block Script
- Location:
- Embroidery Text:
- Embroidery Color:

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:

Please visit www.sunparts.us for spare parts information.