

A quick-folding, lightweight stroller frame with incredibly versatile seating & positioning

EIZ18 (E1232 with seating)

Customer				Ship To			
Account			Name				
Account #		Date		Address			
☐ Quote	☐ Order	PO #		City		State	
Buyer				ZIP		Phone	
Mark For				Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

What's Included

Designed with busy families in mind. The stroller frame quickly and easily folds after removing the seat in one simple step for compact storage.

Multiple seating configurations allows you to personalize your Zippie Voyage to perfectly fit your needs. Each unit includes:

- Foot-Operated Wheel Lock
- Quick-Release Wheels
- 35° Tilt
- 12" Rear Tires & 8" Front Casters
- Angle-Adjustable Stroller Handle
- Forward & Rearward Facing Seating Capability
- Front Caster Suspension & Swivel Locks
- Adjustable Seating Shell
- Under Carriage Storage
- 75 lbs. Weight Capacity
- Canopy with Breathable Window & Mesh Pockets
- (may be limited when some options chosen)

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
Voyage Frame Style			
Frame Seating Options			\$
☐ 005BS7	E1232	Voyage for Advanced Seating with Tilt (includes frame and seating)	\$ 2,912.00
☐ 005BS9	E1232	Voyage for Moderate Seating with Tilt (includes frame and seating)	\$ 2,912.00
Booster Base Seating Options			\$
☐ 005BS10	E1399	Voyage Booster Base with Tilt and Advanced Seating	\$ 3,290.00
☐ 005BS11	E1399	Voyage Booster Base with Tilt and Moderate Seating	\$ 3,290.00
☐ 005BS12	E1399	JCM Booster Base with Tilt and Advanced Seating	\$ 4,638.00
☐ 005BS13	E1399	JCM Booster Base with Tilt and Moderate Seating	\$ 4,638.00
Seating Only Options (omit frame)			\$
☐ 005BS15		Advanced Seating with Tilt (omits frame, seating only)	\$ 583.00
☐ 005BS16		Moderate Seating with Tilt (omits frame, seating only)	\$ 583.00
Transit			\$
☐ 005KJ11	K0108	Transit Option	\$ 374.00



For Moderate Seating, See Moderate Seating Section

For Advanced Seating, See Advanced Seating Section

✓ OPTION	DESCRIPTION						RETAIL
Moderate Seating							
Choose Moderate Seating options, prices reflect ordering base which includes seating.							
Moderate Seating Options				\$		\$	
☐ 005KM01	Moderate Seating		NC	☐ 005KM02	E1014	Recline Feature (required)	\$ 408.00
Grab Bar				\$		\$	
☐ 005KR01	K0108	Grab Bar (optional)				\$	169.00
Upholstery Color							NC
Includes inner covers and crossover canopy							
☐ 005KN01	Blue	☐ 005KN02	Red	☐ 005KN03	Beige	☐ 005KN04	Grey
Moderate Seat Only - Positioning and Supports							
Head Supports (optional)				\$		\$	
☐ 005VJ01	E0955	Small Head Supports	\$ 175.00	☐ 005VJ02	E0955	Large Head Supports	\$ 175.00
Medial Supports (optional)				\$		\$	
☐ 005KQ01	E0957	Tall Medial Support	\$ 158.00	☐ 005KQ02	E0957	Contour Medial Support	\$ 164.00
Laterals (optional)				\$		\$	
☐ 005KO01	E0956	Small Laterals (pair)	\$ 175.00	☐ 005KO02	E0956	Large Laterals (pair)	\$ 175.00
Hip Guides (optional)				\$		\$	
☐ 005KP01	E0956	Hip Guides (pair)	\$ 205.00	☐ 005KP02	E0956	Infant Insert (hip guides req)	\$ 147.00
JAY Anterior Pelvic Belts							
Standard 2-Point Belt with 1" Push Button				\$		\$	
☐ 005ZZ11	E0978	4.5" Pad (APB5115)	\$ 106.00	☐ 005ZZ12	E0978	5.5" Pad (APB5116)	\$ 106.00
Standard 2-Point Belt with 1.5" Push Button				\$		\$	
☐ 005ZZ15	E0978	7.5" Pad (APB7121)					
V-Contour 4-Point Belt				\$		\$	
☐ 005ZZ24	E0978	Small Pad with 1" Push Button (APB3716)					
☐ 005ZZ23	E0978	Medium Pad with 1.5" Push Button (APB3727)					
JAY Anterior Trunk Supports							
Classic structure not available with early intervention size							
Anterior Trunk Support Style							\$
☐ 005JY01	E0960	Classic Dynamic (ATS8100 Series)				\$	177.00
☐ 005JY02	E0960	Classic Structure (ATS8300 Series)				\$	177.00
☐ 005JY06	E0960	Center Opening Dynamic (ATS8800 Series)				\$	177.00
☐ 005JY05	E0960	Center Opening Structure (ATS9400 Series)				\$	177.00
Anterior Trunk Support Size			\$		\$		\$
☐ 005JZ05	Early Intervention	NC	☐ 005JZ06	Extra Small	NC	☐ 005JZ01	Small NC
Moderate Seating - Footrests							
☐ 005F329	K0040	Angle Adjustable Footrest (required)				\$	181.00



Grab Bar



Head Supports



Medial Support



Laterals



Hip Guides



Infant Insert

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

Advanced Seating

Choose Advanced Seating options, prices reflect ordering base which includes seating.

Planar and Contoured Seat

Planar and contoured seat cushions include:

- Contoured Medium Density Foam
- 6 - 13" Hip Width
- Built-In Depth Growth
- Removable Neoprene Cover
- Optional 3" Lateral Hip Blocks (Pr)
- Optional Slight Medial Thigh Contour
- Optional Fixed Ischial Ledge




Options			\$
☐ 005CU28		Planar Seat without Hip Blocks	NC
☐ 005CU41	E2292	Planar Seat with Hip Blocks (pair)	\$ 309.00
☐ 005CU42	E2294	Contoured Wedge Seat with Hip Blocks (pair)	\$ 309.00
☐ 005CU43	E2294	Contoured Anti-thrust (ledge) Seat	\$ 309.00
☐ 005CU2		Omit Cushion	NC
Depth (based on 1" t-back)		\$	\$
☐ 005SD78	Short (5 - 8")	NC	☐ 005SD79 Medium (7 - 10") NC
			☐ 005SD80 Long (9 - 13") NC

Adjustable Contour Seat

Adjustable contour seat includes:

- Choice of Foam Density
- Adjustable Pre-ischial Ledge
- Adjustable 1" Tall Foam Wedges for:
 - Lateral Hip and Thigh Support
 - Medial Knee and Thigh Support
- 3" Seat Depth Growth Kit
- Removable Neoprene Cover



Options			\$
☐ 005ST65	E2294	Adjustable Contour Seat (select foam and depth below)	\$ 372.00
☐ 005CU2		Omit Cushion	NC
Foam			NC
☐ 005FM1	Soft Foam	☐ 005FM2 Medium / Soft Foam	☐ 005FM3 Medium Foam
		☐ 005FM8	Omit Foam
Depth (based on 1" t-back)		\$	\$
☐ 005SD78	Short (5 - 8")	NC	☐ 005SD79 Medium (7 - 10") NC
			☐ 005SD80 Long (9 - 13") 4" growth NC

Advanced Seating - Backrest

Backrest Shell

			\$
☐ 005BT111		Short 7" to 13" Height (select desired backrest height setting below, 10" is default)	NC
		☐ 7" ☐ 8" ☐ 9" ☐ 10" ☐ 11" ☐ 12" ☐ 13"	
☐ 005BT112		Tall 10" to 17" Height (select desired backrest height setting below, 13" is default)	NC
		☐ 10" ☐ 11" ☐ 12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐ 17"	

Backrest Option

☐ 005KJ9	E1014	Recline Feature (90 - 135° range)	\$ 393.00
---------------------------------	-------	-----------------------------------	-----------

Backrest Cushion Options

2.5" "T" Back is used for seat depths as small as 3.5". Backrest has removable foam for additional depth growth without adjusting seat pan.

1" Thick T-Back Cushion Options			\$
☐ 005BT35	E2291	Short 8" to 10"	NC
☐ 005BT28	E2291	Medium 10" to 12"	NC
☐ 005BT133		Omit Back Cushion	NC
2.5" Thick T-Back Cushion Options			\$
☐ 005BT115	E2291	Short 8" to 10"	\$ 79.00
☐ 005BT116	E2291	Medium 10" to 12"	\$ 79.00
☐ 005BT133		Omit Back Cushion	NC
☐ 005BT29	E2291	Tall 12" to 14"	NC
☐ 005BT56	E2291	Extra Tall 14" to 16"	NC
☐ 005BT117	E2291	Tall 12" to 14"	\$ 79.00
☐ 005BT118	E2291	Extra Tall 14" to 16"	\$ 79.00

Extra Tall not available with Short Backrest Shell

✓ OPTION	DESCRIPTION		RETAIL
Advanced Seating - Head Support			
Hardware allows the headrest to be flush mounted to backrest when above shell. When choosing a Headrest option, Hardware Bracket must be chosen, Hardware also available individually. Headrest Ht. Adj. Range from seat pan based on selected Backshell Height. (i.e. Short Ht Ranges 9" - 19.5" & Tall Ht Ranges 13" - 23-1/2")			
☐ 005HE40	E0955	Whitmyer 8" Plush	\$ 253.00
☐ 005HE41	E0955	Whitmyer 10" Plush	\$ 253.00
☐ 005HE42	E0955	Whitmyer Small Contour Cradle	\$ 253.00
☐ 005HE43	E0955	Whitmyer Medium Contour Cradle	\$ 253.00
☐ 005HE29	E0955	Whitmyer Small Adjust-A-Plush	\$ 302.00
☐ 005HQ6	K0108	Whitmyer DFS (forhead strap) (squared for single pad, must order whitmyer headrest)	\$ 134.00
☐ 005HE31	E1028	Hardware Bracket Only (can be chosen without headrest)	\$ 118.00

Advanced Seating - Positioning and Supports

For the supports, Voyage will be set up at factory. Insert dimensions on last page spec. drawing. (Not required, will default to std. settings)

Lateral and Anterior Trunk Supports

Lateral brackets height adjustment range based on Backrest Shell Height (For Example: Short Laterals adj. 3" - 12" and Tall Laterals adj. 3" - 17")

Swing-Away Laterals (adjustable 1 - 8" widths)			\$
☐ 005KE41		Swing-Away Adjustable Bracket - Right (optional pad selection below)	\$ 253.00
	E0956	☐ 005KD15 Large Straight Pad - Right (3 x 6")	\$ 112.00
	E0956	☐ 005KD16 Large Curved Pad - Right (3 x 6")	\$ 112.00
	E0956	☐ 005KD13 Small Straight Pad - Right (2 x 6")	\$ 112.00
	E0956	☐ 005KD14 Small Curved Pad - Right (2 x 6")	\$ 112.00
☐ 005KE42		Swing-Away Adjustable Bracket - Left (optional pad selection below)	\$ 253.00
	E0956	☐ 005KD19 Large Straight Pad - Left (3 x 6")	\$ 112.00
	E0956	☐ 005KD20 Large Curved Pad - Left (3 x 6")	\$ 112.00
	E0956	☐ 005KD17 Small Straight Pad - Left (2 x 6")	\$ 112.00
	E0956	☐ 005KD18 Small Curved Pad - Left (2 x 6")	\$ 112.00
Fixed Laterals (adjustable 1 - 10" widths)			\$
☐ 005KE24		Fixed Angle Adjustable Bracket - Right (optional pad selection below)	NC
	E0956	☐ 005KD15 Large Straight Pad - Right (3 x 6")	\$ 112.00
	E0956	☐ 005KD13 Small Straight Pad - Right (2 x 6")	\$ 112.00
☐ 005KE25		Fixed Angle Adjustable Bracket - Left (optional pad selection below)	NC
	E0956	☐ 005KD19 Large Straight Pad - Left (3 x 6")	\$ 112.00
	E0956	☐ 005KD17 Small Straight Pad - Left (2 x 6")	\$ 112.00

JAY Anterior Trunk Support - Contour			\$
☐ 005A827	E0960	JAY Contour Dynamic Early Intervention Pad	\$ 170.00
☐ 005A828	E0960	JAY Contour Dynamic Extra Small Pad	\$ 170.00
☐ 005A829	E0960	JAY Contour Dynamic Small Pad	\$ 170.00
JAY Anterior Trunk Support - Classic			\$
☐ 005K11	E0960	JAY Classic Dynamic Early Intervention Pad	\$ 177.00
☐ 005K12	E0960	JAY Classic Dynamic Extra Small Pad	\$ 177.00
☐ 005K13	E0960	JAY Classic Dynamic Small Pad	\$ 177.00
☐ 005SS1	E0960	Shoulder Strap Pads (only available with classic trunk support)	\$ 29.00

Anterior Trunk H Strap Support			\$
☐ 005SS2	E0960	XX Small H Strap 5-Point	\$ 185.00

Lower Body Supports

Lateral Hip Support (adjusts from 7.5 - 11" width)			\$
--	--	--	----

Only available with planar seat without hip blocks and adjustable contour seat

☐ 005KC6	E0956	Small - Right (4 x 4", 3/8" thick)	\$ 112.00
☐ 005KC4	E0956	Small - Left (4 x 4", 3/8" thick)	\$ 112.00
☐ 005KC5	E0956	Large - Right (5-3/8 x 4.75", 1/2" thick)	\$ 112.00
☐ 005KC3	E0956	Large - Left (5-3/8 x 4.75", 1/2" thick)	\$ 112.00

✓ OPTION	DESCRIPTION				RETAIL	
Lower Body Supports (continued)						
Lateral Thigh Support (adjusts from 3.5 - 11" width)					\$	
Only available with planar seat without hip blocks and adjustable contour seat						
☐ 005KF6	E0953	Small - Right (4 x 4", 3/8" thick)			\$ 112.00	
☐ 005KF4	E0953	Small - Left (4 x 4", 3/8" thick)			\$ 112.00	
☐ 005KF5	E0953	Large - Right (5-3/8 x 4.75", 1/2" thick)			\$ 112.00	
☐ 005KF3	E0953	Large - Left (5-3/8 x 4.75", 1/2" thick)			\$ 112.00	
Medial Support					\$	
☐ 005KH12	E0957	Small (3"H x 2.5"W x 4.5"L)			\$ 141.00	
☐ 005KH13	E0957	Large (3.75"H x 4"W x 6"L)			\$ 141.00	
☐ 005HE77	E1028	Hardware Medial Support Bracket			\$ 253.00	
JAY Pelvic Positioning Belt					\$	
☐ 005A14	E0978	1" Side Release Padded (APB6121)			\$ 55.00	
☐ 005A15	E0978	1.5" Side Release Padded (APB6222)			\$ 55.00	
☐ 005A263	E0978	1" Push Button Padded (APB5125)			\$ 65.00	
☐ 005A256	E0978	1.5" Push Button Padded (APB7221)			\$ 65.00	
☐ 005A830	E0978	Medial 4-Point Belt with Extra Small Pad			\$ 106.00	
☐ 005A831	E0978	Medial 4-Point Belt with Small Pad			\$ 106.00	
Advanced Seating - Footrest						
Footrest adjustment range is 1.5 - 12". Choose desired preset range - if not selected, long range is default. For the short range, the footplate will be inverted.						
☐ 005F312	K0040	Angle-Adjustable Footrest (required)			\$ 173.00	
☐ 005F313	K0037	Preset in Short Range (1.5 - 4")			NC	
☐ 005F314		Preset in Long Range (4.5 - 12")			NC	
Advanced Seating - Upholstered Colors						
Solid Color Options		\$		\$		
☐ 005CL9	Purple	NC	☐ 005CL13	Orange	NC	
☐ 005CL20	Pink Sorbet	NC	☐ 005CL17	Peacock	NC	
☐ 005CL11	Green	NC	☐ 005CL16	Black	NC	
☐ 005CL12	Navy Blue	NC				
Print Color Options		\$		\$		
Standard trim listed behind each print, if alternative solid color is selected, personalized option is required						
☐ 005CL19	Leopard (standard with black trim)	NC	☐ 005CL15	Red Race Flags (standard with red trim)	NC	
☐ 005CL14	Butterflies (standard with purple trim)	NC				
Canopy Pattern						\$
☐ 005UF01	Cross Over				NC	
Personalized Option						\$
Only available with prints						
☐ 005C001	Personalized Color Selection: Choose and Enter from Solid Color Options				\$ 79.00	
Accessories and Options						
Vent and Battery Trays						\$
☐ 005LF1	E1029	Dovetail Vertical Frame Mounted Hardware for LTV Ventilators (15lb limit, includes 1 utility hook)			\$ 420.00	
☐ 005LF11	E1029	Dovetail Vertical Frame Mounted Hardware (Female) for Trilogy Ventilators (15lb limit, includes 1 utility hook)			\$ 420.00	
☐ 005LF2	E1029	Battery & Vent / Utility Underframe Tray (50lb limit, ideal for U1 or group 22 batteries, includes 1 utility hook)			\$ 459.00	
☐ 005LF3	E1029	Front Mount Utility Platform (50lb limit)			\$ 420.00	
☐ 005LF4	E1029	Rear Add-on Accessory Platform (45lb limit, includes 1 utility hook, must be removed or folded before folding frame)			\$ 196.00	
Upper Extremity Support Tray						\$
Angle-adjustable hardware included. Tray mount adjusts from 1.5 - 7.5" with 1" tall cushion. Cutout is 6")						
☐ 005A91	E0950	Light Grey			\$ 163.00	
☐ 005A92	E0950	Clear Poly			\$ 185.00	

✓ OPTION	DESCRIPTION		RETAIL
Accessories and Options (continued)			
Floor Sitter / Feeder Base			\$
☐ 005A083	E1399	Folding Booster Base (height adjustable)	\$ 1,564.00
☐ 005A082	E1399	JCM Folding Booster Base (height adjustable)	\$ 2,908.00
Shoe Holders and Foot Straps			\$
For moderate seating, footplate does not have mounting slots for shoe holders.			
☐ 005A68	E0952/51	6" Long Shoe Holders (pair)	\$ 134.00
☐ 005A69	E0952/51	8" Long Shoe Holders (pair)	\$ 134.00
☐ 005A79	E0952	16" Long Foot Straps with Covers (pair)	\$ 40.00
☐ 005A80	E0952	18" Long Foot Straps with Covers (pair)	\$ 40.00
Accessories			\$
☐ 005A74	K0105	IV Support Pole	\$ 163.00
☐ 005A713	E2208	O2 Tank Holder	\$ 196.00
☐ 005A714	K0105/E2208	IV / O2 Combination	\$ 274.00
☐ 005LF5		Cup / Bottle / Feeding Holder	\$ 50.00
☐ 005LF7		Utility Hook (i.e. suction canister)	\$ 29.00
☐ 005A462		Rain Cape	\$ 196.00



Front Mount Utility Platform



Cup Holder



Utility Hook



IV Support Pole



Upper Extremity Support Tray



Folding Booster Base



JCM Folding Booster Base



O2 Tank Holder



Battery and Vent/ Utility Underframe Tray




Dovetail Vertical Vent Mount (LTV/ Trilogy)



Rear Add On Accessory Platform

Child Measurements (optional)

Child's Name				1. Top of Head to Seat	
Age		Diagnosis		2. Top of Shoulder to Seat*	
Height			Weight	3. Under Arm to Seat	



4. Trunk Depth
5. Actual Sitting Depth*
6. Foot Length
7. Width of Head
8. Shoulder Width
9. Chest Width
10. Hip Width
11. Left Lower Leg Length*
12. Right Lower Leg Length*

*Critical Dimensions