SUNRISE MEDICAL

QUICKIE®
Live without limits.

ORDER FORM



To see QUICKIE Power Options or Accessories, scan here or visit: https://www.sunrisemedical.com/power-wheelchair-options-accessories


Q500M with SEDEO® PRO™ and JAY® Syncra™ Seating

EIPW30

Customer		Ship To			
Account #		Name			
Date		Address			
PO#		City			
Buyer		State		Zip	
Mark For		Phone			

Client Measurements

a)	Top of head to seat	in
b)	Top of shoulder to seat	in
c)	Underarm to seat	in
d)	Elbow to seat	in
e)	Top of leg to seat	in
f)	Shoulder width	in
g)	Chest width	in
i)	Hip width	in
j)	Actual sitting depth	in
m)	Foot length	in
n)	Knee to foot	in
o)	Seat width	in



Notes

Chair Specifications

HCPC	K0848 / K0856 / K0861
Weights starting at: (with Batteries)	NPC: 300 lbs SPC: 324 lbs MPC: 338 lbs
Power Seating	50° of CG Tilt and Power Recline with 4" mechanical sheer reduction / 12" seat elevator
Seat Type	SEDEO® PRO™ / JAY® Syncra™
Seat Size Range	Width 16" to 22" Depth 14" to 22"
Seat-to-Floor Height Range	16.5" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Mid-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 24", Length 33"
Minimum Turning Radius	18.5" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A
Driving Range	Up to 13.5 miles
Drive Motors	4-pole 6 mph / 10 kph
Drive Wheels	13" w/flat-free inserts (black)
Casters	7" front / 7" rear (black)
Charger	8A Dual Mode (Compact)
Battery Size	Group 34 / Group 22
Color	8 shroud colors and 6 color accents
Ground Clearance	3"

BASE

OPTION	HCPC	DESCRIPTION	RETAIL
Group 3 Standard Duty Mid Wheel Drive Base *End Item Number EIPW30 (PRO™ Seat)*			$
☐ 300PC215	K0848	Quickie Q500 M NPC No Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,541.00
☐ 300PC216	K0856	Quickie Q500 M SPC Single Power Option for SEDEO® PRO™ Seat Included in Base Pricing : Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,541.00
☐ 300PC217	K0861	Quickie Q500 M MPC Tilt/Recline for SEDEO® PRO™ Seat Included in Base Pricing: Q500 M Base, VR-2 90amp Electronics for No Power Seat Functions, Six-Wheel Suspension, 13" Drive Wheels with Dual Density Inserts, 7" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 6mph, Weight Capacity = 300lbs., STFH = 16.5".	$8,541.00
User Weight			$
☐ 300WG23	–	Up to 300 lbs.	NC
Motor Package			$
☐ 300PF99	–	4-pole 10 kph / 6 mph *(Standard for Group 3 K0848, K0856, K0861)*	NC
Tracking			$
☐ 300MI01	K0108	SureTrac® Module & Motor Encoders *(Improves drive tracking)*	$1,628.00
Drive Wheels			$
☐ 300DW69	–	Drive Wheel, 13", Dual Density, Knobby, Black	NC
☐ 300DW58	–	Drive Wheel, 14", Dual Density, Knobby, Black	NC
Caster Wheels			$
☐ 300CW59	–	Caster Wheels 7", Solid, Black	NC
Shroud Color NEW *Base Shroud and Back Shroud are painted to match.*			$
☐ 300C119	–	Matte Black	NC
☐ 300C165	–	Candy Apple Red	NC
☐ 300C166	–	Candy Blue	NC
☐ 300C167	–	Candy Purple	NC
☐ 300C176	–	Green Apple	NC
☐ 300C178	–	Pearl Pink	NC
☐ 300C182	–	Silver	NC
☐ 300C170	–	American Flag	$602.00
Color Accents *Note: Accent colors and shroud colors may not match.*			$
☐ 300EZ20	–	Apple Green	NC
☐ 300EZ22	–	Black	NC
☐ 300EZ21	–	Blue	NC
☐ 300EZ23	–	Orange	NC
☐ 300EZ24	–	Red	NC
☐ 300EZ25	–	Silver	NC
☐ 300EZ19	–	White	NC

BASE

OPTION	HCPC	DESCRIPTION	RETAIL
Service Decals			$
☐ 300D10	–	Army	NC
☐ 300D11	–	Air Force	NC
☐ 300D13	–	Marines	NC
☐ 300D14	–	Navy	NC
☐ 300D15	–	Purple Heart	NC
ISO Transit Tie Downs			$
☐ 300TO50	–	ISO Transit Package	NC
LED Lights			$
☐ 300LI08	–	LED Fender Lights	$566.00
☐ 300LI09	–	LED Fender Lights *(Available with multi-power only)*	NC
☐ 300LI02	–	Omit LED Lights	NC
Seat Heights			$
		Seat Height *Measured from front of seat pan to floor with 0° frame angle.*	
☐ 300HT103	–	16.5"	NC
☐ 300HT55	–	17"	NC
☐ 300HT104	–	17.5"	NC
☐ 300HT56	–	18"	NC
☐ 300HT105	–	18.5"	NC
☐ 300HT57	–	19" *(Only available with Tilt or Seat Elevator.)*	NC
Batteries and Chargers			$
Group 34 Battery Options			
☐ 300AG27	E2359	GP34 Flat Top Pair w/Straps MK Gel	$1,097.00
☐ 300AG29	E2359	GP34 Flat Top Pair w/Straps AGM	$828.00
Group 22 Battery			
☐ 300AG31	E2361	GP22 Flat Top Pair w/Straps	$613.00
Omit Batteries			
☐ 300AG24	–	Battery Harnesses Only *(Omit Batteries)*	NC
Chargers			
☐ 300BJ30	–	10A Dual Mode (Compact)	$110.00
☐ 300BJ19	–	8A Dual Mode (Compact)	NC

POWERED SEATING

	OPTION	HCPC	DESCRIPTION	RETAIL
Power Seating Options				$
Seat Elevator				
☐	300PI07	E2298	Seat Elevator, Powered, 300 mm / 12"	$3,791.00
☐	300PI08	E2298	C-Me® Seat Elevator, Powered, 300 mm / 12" *(C-Me speed up to 3 mph) (N/A with VR-2)*	$3,791.00
Seat Tilt				
☐	300PA48	E1002	Posterior 50° Tilt	$8,458.00
Seat Recline				
☐	300PA37	E1004	150° Powered Recline w/ 4" Shear Reduction	$8,458.00
Seat Tilt/Recline				
☐	300PA41	E1007	Posterior 50° Tilt and 150° Recline Powered, w/ 4" Shear Reduction	$13,709.00
Power Recline Backrest Preset Home Angle				
☐	300QK1	–	90°	NC
☐	300QK2	–	95° (Default Position)	NC
Pre Set Seat Frame Angle *(No Power Seating Options, Increases Seat to Floor Height)*				$
☐	300FA10	–	0° *(Adds 0" to STFH)*	NC
☐	300FA3	–	3° *(Adds 1" to STFH)*	NC
Backrest Frame Type				$
☐	300BF16	–	Tool Angle Adjustable Back, 90°-120°	NC
☐	300BF14	E1225	User Adjustable Manual Semi-Recline Back, 90°-130°, Right *(Adjustment Lever Mounted on Opposite Side of Joystick)*	$682.00
☐	300BF15	E1225	User Adjustable Manual Semi-Recline Back, 90°-130°, Left *(Adjustment Lever Mounted on Opposite Side of Joystick)*	$682.00
Backrest Frame Option				$
☐	300BF13	–	Flip Down, Quick Release	$147.00
Back Cushion Type				$
☐	300BT134	–	SEDEO® Back Cushion *(Choose Seating in Section A)*	NC
☐	300BT138	–	JAY® Syncra™ *(Choose Seating in Section B)*	NC
☐	300BT135	–	JAY® J3® *(Choose Seating in Section C)*	NC
☐	300BT137	–	SEDEO® PRO™ Aftermarket Back Adapter *(Height Adjustable Round Tubes for Traditional Mounting)*	$525.00
☐	300BT133	–	OMIT *(For personalized Mounting. No Headrest Mount)*	NC

POWERED SEATING

	OPTION	HCPC	DESCRIPTION				RETAIL
A - SEDEO® Backrest Dimensions							$
Back Width							
☐	300BN28	–	400mm / 16"				NC
☐	300BN29	–	460mm / 18"				NC
☐	300BN30	–	500mm / 20"				NC
☐	300BN31	–	560mm / 22"				NC
Back Shell Height				**Tool Angle Adjust.**	**150° Recline**		
☐	300B255	–	460mm / 18"	22"-23"	23"-25"	**Installed Back Height**	NC
☐	300B256	–	560mm / 22"	24"-27"	26"-29"		NC
☐	300B257	–	660mm / 26"	28"-30"	30"		NC
Installed Back Height - *(From Pan to Top of Back)*							
☐	300JR32	–	22"				NC
☐	300JR33	–	23"				NC
☐	300JR34	–	24"				NC
☐	300JR35	–	25"				NC
☐	300JR36	–	26"				NC
☐	300JR37	–	27"				NC
☐	300JR38	–	28"				NC
☐	300JR39	–	29"				NC
☐	300JR40	–	30"				NC

B - JAY® Syncra™ Seating

B - JAY® Syncra™ Back P. 21 $

☐	300JP14	E2620	JAY® Syncra™ Back 16-20" Wide with 3" Contour Depth	$1,138.00
☐	300JP15	E2621	JAY® Syncra™ Back 22" Wide with 3" Contour Depth	$1,255.00

B - JAY® Syncra™ Back Width $

Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.

☐	300JQ07	–	16" Wide *(Adjustment range 16 to 17" Wide)**	NC
☐	300JQ08	–	18" Wide *(Adjustment range 18 to 19" Wide)**	NC
☐	300JQ09	–	20" Wide *(Adjustment range 20 to 21" Wide)**	NC
☐	300JQ10	–	22" Wide *(Adjustment range 22 to 23" Wide)*	NC

**Additional +1" of shell width adjustment is achievable. Requires purchase of 1" larger foam and upholstery separately.*

B - JAY® Syncra™ Back Shell Options and Height Ranges $

JAY® Syncra™ Back shells are mounted with a minimim 5.5" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field with-in the ranges shown based on the back frame type.

JAY® Syncra™ Shell Height

☐	300JR42	–	16" Tall Shell *(22" to 24" Back Height)*	NC
☐	300JR43	–	18" Tall Shell *(24" to 26" Back Height)*	NC
☐	300JR44	–	20" Tall Shell *(26" to 28" Back Height)*	NC
☐	300JR45	–	22" Tall Shell *(28" to 30" Back Height)*	NC

POWERED SEATING

	OPTION	HCPC	DESCRIPTION	RETAIL
JAY® Syncra™ Back Installed Back Height (From Pan to Top of Back)				
☐	300JR50	–	22" *(16" Shell Height)*	NC
☐	300JR51	–	23" *(16" Shell Height)*	NC
☐	300JR52	–	24" *(16" Shell Height)*	NC
☐	300JR53	–	24" *(18" Shell Height)*	NC
☐	300JR54	–	25" *(18" Shell Height)*	NC
☐	300JR55	–	26" *(18" Shell Height)*	NC
☐	300JR56	–	26" *(20" Shell Height)*	NC
☐	300JR57	–	27" *(20" Shell Height)*	NC
☐	300JR58	–	28" *(20" Shell Height)*	NC
☐	300JR67	–	28" *(22" Shell Height)*	NC
☐	300JR68	–	29" *(22" Shell Height)*	NC
☐	300JR69	–	30" *(22" Shell Height)*	NC
B - JAY® Syncra™ Back Upholstery				$
☐	300BU96	–	Syncra™ Premium Leatherette	NC
☐	300BU97	–	Syncra™ Stretch	NC
B - JAY® Syncra™ Fixed Lateral Support Accessory				$
JAY® Fixed Lateral Support Accessory Hardware *(JAY® Fixed Lateral Support is height adjustable)*				
☐	300JW33	E0956	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Right	$147.00
☐	300JW34	E0956	Fixed Lateral Support Accessory (L x H) 90 mm x 140 mm / 3.5" x 8" - Left	$147.00
C - JAY® J3® Back				
C - JAY® J3® Back - Contour				$
☐	300JP01	E2613	J3® PA Back for 16"-20" Wide	$625.00
☐	300JP02	E2614	J3® PA Back for 22" Wide	$928.00
☐	300JP11	E2615	J3® PL Back for 16"-20" Wide	$729.00
☐	300JP13	E2620	J3® PD Back for 16"-20" Wide *(18" and 21" Shell Height Only)*	$875.00
C - JAY® J3® Back Width				$
Back can be ordered a different width than seat pan. Armrest will be adjusted to back width ordered.				
☐	300JQ02	–	16"	NC
☐	300JQ03	–	18"	NC
☐	300JQ04	–	20"	NC
☐	300JQ05	–	22"	NC

POWERED SEATING

	OPTION	HCPC	DESCRIPTION				RETAIL
C - JAY® J3® Back Shell Options and Height Ranges							$
			JAY® J3® shells are mounted with a minimim 4" gap between the bottom of the shell to the seat pan to accommodate the mechanical shear. Back height can be adjusted in the field within the ranges shown.				
JAY® J3® Shell Height				**Tool Angle Adjust.**	**150° Recline**		
☐	300JR41	–	15" tall shell *(Upper Thoracic Med) for 22" Wide*	22"-23"	22"-23"	Installed Back Height	NC
☐	300JR14	–	16.5" tall shell *(Upper Thoracic Short)*	21"-23"	22"-24"		NC
☐	300JR15	–	18" tall shell *(Upper Thoracic Med)*	22"-24"	23"-25"		NC
☐	300JR16	–	19.5" tall shell *(Upper Thoracic Tall)*	24"-26"	24"-26"		NC
☐	300JR18	–	18" tall shell *(Shoulder Height) for 22" wide*	24"	24"		NC
☐	300JR17	–	21" tall shell *(Shoulder Height)*	24"-26"	25"-28"		NC
☐	300JR27	–	22.5" tall shell *(Shoulder Height)*	26"-28"	27"-30"		NC
☐	300JR28	–	24" tall shell *(Shoulder Height)*	27"-29"	28"-30"		NC
JAY® J3® Installed Back Height *(From Pan to Top of Back)*							
☐	300JR29	–	21"				NC
☐	300JR19	–	22"				NC
☐	300JR20	–	23"				NC
☐	300JR21	–	24"				NC
☐	300JR22	–	25"				NC
☐	300JR23	–	26"				NC
☐	300JR24	–	27"				NC
☐	300JR25	–	28"				NC
☐	300JR30	–	29"				NC
☐	300JR31	–	30"				NC
C - JAY® Anterior Thoracic Support							$
			Classic-Structure Pull Down Harness (4 point strap, only available with JAY® backrest)				
☐	300K10	E0960	Medium *(ATS8302)*				$177.00
☐	300K11	E0960	Large *(ATS8303)*				$177.00

POWERED SEATING

	OPTION	HCPC	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad, Upholstery and Mount)				$
JAY® Lateral Thoracic Support Pads				
☐	300JW35	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	$147.00
☐	300JW38	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	$147.00
☐	300JW36	E0956	Lateral Thoracic (L x H) 100 mm x 152 mm / 4" x 6" Pad - Right	$147.00
☐	300JW39	E0956	Lateral Thoracic (L x H) 100 mm x 52 mm / 4" x 6" Pad - Left	$147.00
☐	300JW37	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	$147.00
☐	300JW40	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	$147.00
JAY® Lateral Thoracic Support Pad Upholstery				
☐	300NU01	–	Syncra™ Premium Leatherette	NC
☐	300NU02	–	Syncra™ Stretch	NC
JAY® Lateral Thoracic Swing-Away Mount				
☐	300JW41	E1034	Lateral Thoracic Swing Away Bracket - Right	$292.00
☐	300JW42	E1034	Lateral Thoracic Swing Away Bracket - Left	$292.00
Headrests *(Requires Pad, Upholstery, & Mount)*				$
Headrest Pad				
☐	300HE36	E0955	SEDEO® Pad Small with Cover (LxH) 240 x 150 mm / 10" x 6" *Can only be used with Multi Position Headrest Mount with quick-release and SEDEO® upholstery*	$188.00
☐	300HE38	E0955	Whitmyer® Plush Headrest Pad 203 mm / 8"	$239.00
☐	300HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	$239.00
☐	300HE44	E0955	Whitmyer® Plush Headrest Pad 356 mm / 14"	$250.00
☐	300HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium *(Adjustable 7" to 16" Wide x 4" Tall)*	$367.00
☐	300HE48	–	OMIT Headrest Pad	NC
Headrest Pad Upholstery				
☐	300NU10	–	SEDEO® Upholstery	NC
☐	300NU04	–	Syncra™ Premium Leatherette	NC
☐	300NU11	–	Whitmyer® Lycra / Stretch	NC
Headrest Mount				
☐	300MH32	E1033	SEDEO® Multi-Position Headrest Mount with Quick Release *(Only available with SEDEO® Pad)*	$234.00
☐	300MH33	E1033	Whitmyer® Linx Headrest Mount with Quick Release	$192.00
☐	300MH34	E1033	Whitmyer® Linx-2 Headrest Mount with Quick Release	$234.00
☐	300MH35	E1033	Whitmyer® AXYS Headrest Mount with Quick Release	$123.00
☐	300MH17	–	OMIT Headrest Mount	NC

POWERED SEATING

	OPTION	HCPC	DESCRIPTION	RETAIL
Push Handles				$
☐	300PH27	–	Height Adjustable Push Handles	$350.00
Seat Frame Rails				$
☐	300SN14	–	Seat Rail Short, for Effective Seat Depth 14" - 18"	NC
☐	300SN15	–	Seat Rail Long, for Effective Seat Depth 17" - 22"	NC
Seat Cushion Type				$
☐	300ST81	–	JAY® Syncra™ (See Section D)	NC
☐	300ST80	–	JAY® (See Section E)	NC
☐	300ST75	–	Omit	NC
Seat Size				$
Seat Pan Width				
☐	300R53	–	16"	NC
☐	300R54	–	17"	NC
☐	300R55	–	18"	NC
☐	300R56	–	19"	NC
☐	300R57	–	20"	NC
☐	300R58	–	21" *(Not available with JAY® Syncra™ Cushions)*	NC
☐	300R59	–	22"	NC
Seat Pan Depth ***(Front of seat pan to front of back cushion – effective depth)***				
☐	300RF52	–	14"	NC
☐	300RF53	–	15"	NC
☐	300RF51	–	16"	NC
☐	300RF45	–	17"	NC
☐	300RF46	–	18"	NC
☐	300RF47	–	19"	NC
☐	300RF48	–	20"	NC
☐	300RF49	–	21" *(Not available with JAY® Syncra™ Cushions)*	NC
☐	300RF50	–	22"	NC

POWERED SEATING

	OPTION	HCPC	DESCRIPTION	RETAIL
D - JAY® Syncra™ Cushion P. 21 (JAY Your Way modifications available through Customer Service (800) 333-4000)				$
☐	300CU89	E2622	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology *(16" to 20" wide)*	$630.00
☐	300CU90	E2623	JAY® Syncra™ Cushion featuring JAY Flow® Fluid with Cryo™ Technology *(22" wide)*	$729.00
JAY® Syncra™ Cushion Cover Upholstery				
☐	300YR73	–	Syncra™ Premium Leatherette	NC
☐	300YR74	–	Syncra™ Stretch	NC
E - JAY® Seat Cushion				$
☐	300CU59	E2607	JAY® Ion® *(Sizes 16" to 21" Wide)*	$460.00
☐	300CU60	E2608	JAY® Ion® *(Size 22" Wide)*	$578.00
☐	300CU61	E2607	JAY® Union® *(Sizes 16" to 21" Wide)*	$520.00
☐	300CU62	E2608	JAY® Union® *(Size 22" Wide)*	$635.00
☐	300CU64	E2622	JAY® Fusion® *(Sizes 16" to 21" Wide)*	$514.00
☐	300CU65	E2623	JAY® Fusion® *(Size 22" Wide)*	$613.00

POSITIONING SUPPORTS & ACCESSORIES

	OPTION	HCPC	DESCRIPTION	RETAIL
Knee and Thigh Swing-Away Support Pads				$
Lateral Knee Support Pads and Hardware *(must pick pad and mount)*				
☐	300XY05	E0953	Lateral Knee Support Pad, Right *(4" round)*	$205.00
☐	300XY06	E0953	Lateral Knee Support Pad, Left *(4" round)*	$205.00
☐	300XX08	E1028	Lateral Swing-Away Bracket Knee - Right	$268.00
☐	300XX09	E1028	Lateral Swing-Away Bracket Knee - Left	$268.00
Lateral Thigh Support Pads and Hardware *(must pick pad and mount)*				
☐	300KF22	E0953	Lateral Thigh Support Pad Short - Right *(4.8" H x 4.8 " L / 120mm x 120mm)*	$147.00
☐	300KF23	E0953	Lateral Thigh Support Pad Short - Left *(4.8" H x 4.8 " L / 120mm x 120mm)*	$147.00
☐	300KF24	E0953	Lateral Thigh Support Pad Long - Right *(4.8" H x 7" L / 120mm x 120mm)*	$147.00
☐	300KF25	E0953	Lateral Thigh Support Pad,Long - Left *(4.8" H x 7" L / 120mm x 120mm)*	$147.00
☐	300KF26	E1028	Lateral Thigh Removable Mount - Right	$292.00
☐	300KF27	E1028	Lateral Thigh Removable Mount - Left	$292.00
JAY® Lateral Thigh Supports *(Must pick pad, upholstery, and mount)*				$
JAY® Lateral Thigh Support Pads				
☐	300KF16	E0953	Lateral Thigh Support Pad Standard, Right *(8.25" x 3.25")*	$147.00
☐	300KF17	E0953	Lateral Thigh Support Pad Standard, Left *(8.25" x 3.25")*	$147.00
☐	300KF18	E0953	Lateral Thigh Support Pad Large, Right *(12" x 3.25")*	$147.00
☐	300KF19	E0953	Lateral Thigh Support Pad Large, Left *(12" x 3.25")*	$147.00
JAY® Lateral Thigh Support Pad Upholstery				
☐	300NU07	–	Syncra™ Premium Leatherette	NC
☐	300NU08	–	Syncra™ Stretch	NC
JAY® Lateral Thigh Removable Mount				
☐	300KF20	E1028	Lateral Thigh Removable Mount - Right	$292.00
☐	300KF21	E1028	Lateral Thigh Removable Mount - Left	$292.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Armrest Structure			$
☐ 300EG41	–	Reclining, Flip Back, Short *(Best used with shorter seat depths) (Height adjustable 9" to 15")*	NC
☐ 300EG36	–	Reclining, Flip Back with 2-point post locking arm. *(Height adjustable 11" to 15" x 0.625")*	$292.00
☐ 300EG52	–	Reclining, Flip Back, Dual Link Armrest *(Height adjustable 8" to 15" including length and lateral adjustment)*	NC
☐ 300EG32	–	Single Post, Flip Back *(N/A with thigh support) (Height adjustable 10" to 15")*	NC
Arm Pads			$
SEDEO® Armrest Pad			
☐ 300AB43	–	300 mm / 12", Integral Foam, Right	NC
☐ 300AE43	–	300 mm / 12", Integral Foam, Left	NC
☐ 300AB44	–	400 mm / 16", Integral Foam, Right	NC
☐ 300AE44	–	400 mm / 16", Integral Foam, Left	NC
Comfort Arm Pad *(Foam)*			
☐ 300AB46	–	2" Wide Right	$88.00
☐ 300AE46	–	2" Wide Left	$88.00
Comfort Arm Pad *(Gel)*			
☐ 300AB48	–	2" Wide Right	$175.00
☐ 300AE48	–	2" Wide Left	$175.00
Ottobock® Arm Troughs			
☐ 300AE06	–	Ottobock Arm Trough, Med (12") Left	$136.00
☐ 300AB06	–	Ottobock Arm Trough, Med (12") Right	$136.00
☐ 300AE07	–	Ottobock Arm Trough, Large (14") Left	$136.00
☐ 300AB07	–	Ottobock Arm Trough, Large (14") Right	$136.00
☐ 300OL01	–	Ottobock Hand Pad, Med (7") Left	$59.00
☐ 300OR01	–	Ottobock Hand Pad, Med (7") Right	$59.00
☐ 300OL02	–	Ottobock Hand Pad, Large (9") Left	$59.00
☐ 300OR02	–	Ottobock Hand Pad, Large (9") Right	$59.00
☐ 300AE78	–	Ottobock One Piece Arm Trough (18"), Left	$190.00
☐ 300AB78	–	Ottobock One Piece Arm Trough (18"), Right	$190.00
Elbow Support			$
☐ 3003A02	–	Elbow Support Pad, RH 3" x 10" *(Includes side support)*	$179.00
☐ 3003A03	–	Elbow Support Pad, LH 3" x 10" *(Includes side support)*	$179.00
☐ 3003A06	–	Elbow Support Pad, RH 3" x 3"	$149.00
☐ 3003A07	–	Elbow Support Pad, LH 3" x 3"	$149.00
☐ 3003A04	E1028	Swing-Away Elbow Support Hardware, RH	$245.00
☐ 3003A05	E1028	Swing-Away Elbow Support Hardware, LH	$245.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	HCPC	DESCRIPTION	RETAIL
Side Guards *(Pair)*			**$**
☐ 300SG44	E0953	Side Guards Integral Foam *(11"L x 5.5"H) (N/A with Reclining Armrests)*	$239.00
Fixed Centermount Legrest			**$**
☐ 300H183	–	SEDEO® Seating Fixed Centermount Legrest *(3" to 17.5" Leg Length)*	$292.00
☐ 300H184	–	Fixed Centermount Legrest, Short *(Includes Calf Pads and Single Footplate) (11.5" to 14.5" Leg Length)*	NC
☐ 300H185	–	Fixed Centermount Legrest, Long *(Includes Calf Pads and Single Footplate) (13.5" to 18.5" Leg Length)*	NC
Power Centermount Legrest			**$**
☐ 300H175	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension *(Mechanical Articulation) (5" to 17.5" Leg Length)*	$2,706.00
☐ 300H109	E1012	Power Centermount Legrest, 90° - 0° Angle, 7.5" Extension *(Powered Articulation) (3" to 17.5" Leg Length)*	$3,791.00
Centermount Legrest Accessories			**$**
Footplates			
☐ 300F127	K0040	Dual, Height, Width and Angle Adjustable, Standard *(14"W x 10"D)*	$197.00
☐ 300F128	K0040	Dual, Height, Width and Angle Adjustable, Large *(16"W x 12"D)*	$197.00
Footplate Accessories			
☐ 300DI68	K0108	Extra Large Centermount Footplate *(Single Footplate Overlay 16"W x 12"D)*	$199.00
☐ 300DI69	K0108	Footplate Connector Link *(Attaches footplates together)*	$97.00
☐ 300DI57	K0108	Gel Footplate *(Pair, LH and RH)*	$403.00
Calf Pad Options			
☐ 300DI34	–	Calf Pads, Pair, Standard Integral Foam	NC
☐ 300DI35	–	Calf Pads, Pair, Padded *(N/A w/ Std footplate)*	$408.00
SEDEO® PRO™ Swing-Away Legrest *(Pair) (Footplates, calf supports & options must be selected separately.)*			**$**
☐ 300H177	–	Legrests, Swing-Away, 70°, SEDEO® PRO™ *(14" to 21" Leg Length)*	NC
☐ 300H178	E0990	Legrests, Manual Elevating, Swing-Away, SEDEO® PRO™ *(Includes Padded Calf Pads) (14.5" to 21.5" Leg Length)*	$408.00
☐ 300H179	E1010	Legrests, Power Elevating, Swing-Away, SEDEO® PRO™ *(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)*	$2,508.00
☐ 300DI76	–	Omit Legrests	NC
Swing-Away Legrest Footplates			
☐ 300F369	K0040	Angle Adjustable Flip-Up Footplates	$197.00
Swing-Away Legrest Calf Supports			
☐ 300DI42	–	Calf Support, Depth Adjustable, Padded *(6.5" to 8" Wide, depending on chair width, x 9.5" Tall)*	$460.00
Swing-Away Legrest Options			
☐ 300DI44	E0951	Heel Loops	$58.00
☐ 300KG4	–	Knee Pads, Integral Foam (2.5" Wide x 4.5" Tall)	$234.00
SEDEO® Lower Extremity Swing-Away Support			
☐ 300VD01	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Right	$460.00
☐ 300VD02	E1020	Lower Extremity Support, Padded, Swing-Away, SEDEO®, Left	$460.00

POSITIONING SUPPORTS & ACCESSORIES

	OPTION	HCPC	DESCRIPTION	RETAIL
Pelvic Positioning Belts				
☐	300A837	–	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
☐	300A839	–	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
☐	300A840	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	$88.00
☐	300A845	K0108	Bodypoint – Evoflex Pelvic Stabilizer Belt	$203.00
☐	300A846	K0108	Bodypoint – Evoflex Pelvic Stabilizer Belt with Extension	$224.00
Accessories				$
Ventilator Option				
☐	300VT12	E1029	Vent Tray Fixed *(Pulmonetic LTV Models)*	$542.00
O2 Tank Option				
☐	300A718	E2208	O2 Tank Holder Mounted to Backrest *(D and E Styles)*	$379.00
Wheelchair Tray Table				
☐	3007A09	E0950	Tray, Swing-Away	$408.00
Backpack				
☐	300A1	–	Black Utility	$63.00
Medial Thigh Support				
☐	3008A80	E0957	Medial Thigh Support	$197.00
☐	3008A81	E1028	Pommel Swing-Away Hardware	$292.00
Misc. Accessories				
☐	3008A59	–	Transfer Bar *(Pair)*	$408.00
☐	3008A52	–	Ram Mount Phone	$100.00
☐	3008A53	–	Ram Mount Tablet	$408.00
☐	3008A54	–	Ram Mount Cup Holder	$100.00
☐	3008A58	–	Rear View Mirror	$100.00
☐	300A611	–	Crutch Holder	$147.00
☐	3008A73	–	Bag For Armrest	$100.00
☐	3008A75	–	Reflectors *(Backrest shroud mounted)*	$100.00
☐	3009A10	–	Utility Hooks	$60.00

JOYSTICK / ELECTRONICS

VR-2 Joystick			$
☐	300J127	–	VR-2 Joystick	NC

R-net Joystick *(Remote)*			$
☐	300J121	–	Standard LED Joystick *(w/o Screen Display)*	NC
☐	300J116	–	Color Joystick *(2.3" Color LCD Screen Display, CJSM1)*	NC
☐	300J118	–	High Resolution Color Joystick, (CJSM2+BT) *3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control*	$1,230.00
☐	300J119	–	Omit Remote Control *(Required with arm mounted specialty controls)*	NC

Joystick Mounting			$
☐	300JM34	–	Fixed on Armrest	NC
☐	300JM35	E1032	Parallel Swing-Away	$345.00
☐	300JM36	E1032	Link-It Mid Line Mount on Armrest *(Swing-Away)*	$991.00
☐	300JM49	–	Attendant *(Available with LED or CJSM1 only. Mounted to Push Handle)*	NC

Joystick Mount (Location)			$
☐	300HM16	–	Right Hand	NC
☐	300HM17	–	Left Hand	NC

Display Language Selection R-net			$
☐	300L1	–	English	NC
☐	300L3	–	French	NC
☐	300L7	–	Spanish	NC

Joystick Handle (Optional)			$
☐	300JH13	E2323	Ball Handle	$147.00
☐	300JH15	E2323	T-Handle	$147.00
☐	300JH24	E2323	Bodypoint – Rubber Dome Joystick Handle	$191.00
☐	300JH25	E2323	Bodypoint – U-Shaped Joystick Handle *(4 in)*	$191.00
☐	300JH26	E2323	Bodypoint – Grooved Mushroom Joystick Handle	$191.00

Electronic Programmers			$
☐	300PR13	–	R-net PC Wireless Programmer	$613.00

R-net Power Modules			$
☐	300ER12	–	Power Module, VR-2 90A, Non-Expandable	NC
☐	300ER22	–	Power Module, R-net, 90A, Non-Expandable *(Only available with single power)*	NC
☐	300ER13	–	Power Module, R-net, 90A, Expandable *(Required for Specialty Controls or 2+ Actuators)*	$1,335.00
	300EX11	E2377	Expandable 90A Controller - R-net	$706.00
	300EX12	E2313	Harness for 90A Expandable Controller	$629.00

JOYSTICK / ELECTRONICS

OPTION	HCPC	DESCRIPTION	RETAIL
Seating Control			$
Powered Seating Control			
☐ 300EO40	–	Omit Thru-Drive *(Requires Switch Control)*	NC
☐ 300EO37	E2310	Thru-drive Control for 1 Actuator	$2,240.00
☐ 300EO38	E2311	Thru-drive Control for 2+ Actuators	$4,055.00
Switch Control *(Available with Omit Thru-Drive only) (1 Actuator for Tilt or Elevate only)*			
☐ 300EO33	–	Dual Toggle for 1 Actuator	NC
☐ 300EO34	–	2.5" Single Button Switch for 1 Actuator	NC
☐ 300EO41	–	Ctrl+5 Switch Box with Voice Activation for 2 Actuators or more	NC
Additional Power Seating Control *(Available with Thru-Drive only)*			
☐ 300QG11	–	Ctrl+5 Switch Box with Voice Activation	$1,459.00
☐ 300QG12	–	Ctrl+5 Switch Box with Voice Activation *(Multi-power only)*	NC
Dual Toggle Mounting Location			
☐ 300TN1	–	Toggle Mount, Right	NC
☐ 300TN2	–	Toggle Mount, Left	NC
Ctrl+5 Mounting Location *(Will swing away if mounted on the same side as the joystick)*			
☐ 300EY19	–	Right Armrest	NC
☐ 300EY20	–	Left Armrest	NC
☐ 300EY21	–	Back Cane, Right	$175.00
☐ 300EY24	–	Back Cane, Left	$175.00
Ctrl+5 Switch Type			
☐ 300EY22	–	Push Buttons	NC
☐ 300EY23	K0108	Toggle Switches	$350.00
Attendant Controls			$
☐ 3000A05	E2331	Dual Control Module for Attendant	$1,517.00
☐ 3000A06	E2331	R-net Color LCD (CJSM1), Control Attendant Only	$1,517.00
Attendant Mounting Location			$
☐ 300DN5	–	Right Back Cane	NC
☐ 300DN6	–	Left Back Cane	NC

JOYSTICK / ELECTRONICS

	OPTION	HCPC	DESCRIPTION	RETAIL
R-net Programming				$
	Assignable Buttons			
	Assignable Button feature allows for personalization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.			
☐	300PM12	–	Include Assignable Buttons	NC
	Power Centermount Footplate Extension Option			
☐	300IE03	K0108	Include 7.5" Powered Footplate Extension for Transfer Assist *(Electronics programmed to allow footplates to be extended with or without legrest being raised.) Only Available with Powered Articulation*	$1,459.00
R-net Omni2 Options	*For Expandable Controls, choose all that apply.*			$
☐	300ON14	–	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth® Choose Mount Below	$2,072.00
☐	300ON15	–	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared and Bluetooth® Bundled Omni2 Available w/ Switch-It Input Devices	NC
Omni2 Mount	*If used w/ Standard Joystick, it must be mounted opposite side.*			$
☐	300DS14	E1032	Link-It with Parallel Swing-Away Omni2 Mount w/Rear Mounted Interface, Right	$578.00
☐	300DS15	E1032	Link-It with Parallel Swing-Away Omni2 Mount w/Rear Mounted Interface, Left	$578.00

SPECIALTY INPUT DEVICES

	OPTION	HCPC	DESCRIPTION	RETAIL
Hand Control Packages *Requires Standard Joystick Omit / Only one hand control device can be ordered.*				$
Proportional Input Devices				
☐	300EU14	E2312	Switch-It MicroGuide	$3,844.00
☐	300EU17	E2328	Switch-It VersaGuide	$2,118.00
☐	300EU18	E2328	Switch-It VersaGuide EZ	$2,502.00
Mounting Position *(Required)*				
☐	300EM04	–	Link-It to the Backrest *(Not Available w/ Omit Back Cushion Type)*	$1,109.00
☐	300EM05	–	Link-It to the Fixed Joystick Mount	$345.00
☐	300EM06	E1032	Link-It to the Parallel Swing-Away Joystick Mount	$578.00
☐	300EM07	E1032	Link-It to the Midline Swing-Away	$1,225.00
☐	300EM08	–	Omit Mount	NC
Mounting Platforms *(Optional)* *Not available w/ Link-It to Backrest*				
☐	300MG16	E0950	SunDisk Hand Tray *(Only available w/ MicroSeries Joysticks)*	$345.00
☐	300MG15	E0950	Bullet Hand Tray	$345.00
Mounting Location *(Required)*				
☐	300IM09	–	Right	NC
☐	300IM10	–	Left	NC
Switch Options				
☐	–	–	See Switch Selections *(Maximum 2)*	–
Head Control Packages *(Select headrest mount on next page)*				$
Only one head control device can be ordered.				
Proportional Input Devices				
☐	300EU25	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads	$6,658.00
☐	300EU26	E2328	Switch-It DualPro® Proportional Head Array - Spot Pads	$6,658.00
☐	300EU51	E2328	Switch-It Vigo Head Drive	$5,648.00
Switched Input Devices				
☐	300EU28	E2330	Switch-It 3-Switch Head Array - Short Curved Wings	$5,127.00
☐	300EU30	E2330	Switch-It 3-Switch Head Array - Long Curved Wings	$5,127.00
Switch Options				
☐	–	–	See Switch Selections *(Maximum 2)*	–

SPECIALTY INPUT DEVICES

	OPTION	HCPC	DESCRIPTION	RETAIL
Sip 'N Puff				$
☐	300ID73	E2330	Switch-It Sip 'N Puff Combo Head Array with Link-It Hardware *(Requires headrest mount)*	$5,827.00
☐	300ID74	K0108	Switch-It OptiStop with Sip 'N Puff	$2,270.00
Headrest Mounts *(Applies to Head/Sip 'N Puff control packages)*				$
☐	300IJ13	E1033	Whitmyer® Linx Headrest Mount with Quick Release	$192.00
☐	300IJ14	E1033	Whitmyer® Linx-2 Headrest Mount with Quick Release	$234.00
☐	300IJ15	E1033	Whitmyer® AXYS Headrest Mount with Quick Release	$123.00
Chin Control Packages				$
Proportional Input Devices				
☐	300EU34	E2312	Switch-It MicroGuide *(Switches mounted to the headrest only.)*	$3,844.00
☐	300EU36	E2328	Switch-It VersaGuide	$2,118.00
☐	300EU37	E2328	Switch-It VersaGuide EZ	$2,502.00
Mounting Position				
☐	300EM09	E1032	Link-It To Chin Harness, Small (150mm / 6")	$1,459.00
☐	300EM10	E1032	Link-It To Chin Harness, Large (220mm / 8")	$1,459.00
☐	300EM11	E1032	Link-It To Headrest, Swing-Away	$1,459.00
☐	300EM16	–	Omit Mount	NC
Input Device Mounting Location				
☐	300IM13	–	Right	NC
☐	300IM14	–	Left	NC
Switch Options				
☐	–	–	See Switch Selections	–
Safety Devices				$
☐	300SA06	–	Switch-It Remote Stop Switch - XLR Connector	$1,091.00
Switched Driving Packages				$
☐	300SW44	E2328	Switch-It™ Cool Cube with DB9 connector for proximity/proportional switch driving. *(Pro Spot proportional switches are required if proportional driving is desired.)* *(For non-proportional driving, single switches can be used.)*	$4,048.00

SPECIALTY INPUT DEVICES

OPTION	HCPC	DESCRIPTION	RETAIL

Switches

Single switches can be programmed to activate any function of the wheelchair including: On/Off, Mode, Profile, Speeds, Seating, Memory Positions, Horn, Lights & Indicators, IR functions, Bluetooth mouse clicks, Environmental Controls when plugged into assignable port.
Switches can connect via: Joysticks (Standard, Micro series, Versa series) Omni2, head controls
Switch ordering instructions:

1. Select switch type and quantity in Section "A" using colums below. (Max 6 switches. See limitations below.)
2. Select switch mounting position from Section "B" using corresponding switch column.
3. Select switch mounting location from Section "C" using corresponding switch column.

Section A – Switch Selections — **Switch / Quantity**

Single Function Switches - Mechanical Switches	Switch 1	Switch 2	Switch 3	Switch 4	Switch 5	Switch 6	$
2.5" Single Switch - Black	☐	☐	☐	☐	☐	☐	$227.00
2.5" Single Switch - Blue	☐	☐	☐	☐	☐	☐	$227.00
2.5" Single Switch - Green	☐	☐	☐	☐	☐	☐	$227.00
2.5" Single Switch - Red	☐	☐	☐	☐	☐	☐	$227.00
2.5" Single Switch - Yellow	☐	☐	☐	☐	☐	☐	$227.00
Egg Switch - Black	☐	☐	☐	☐	☐	☐	$129.00
Egg Switch - Blue	☐	☐	☐	☐	☐	☐	$129.00
Egg Switch - Green	☐	☐	☐	☐	☐	☐	$129.00
Egg Switch - Red	☐	☐	☐	☐	☐	☐	$129.00
Egg Switch - Yellow	☐	☐	☐	☐	☐	☐	$129.00
MicroLite *(Includes color labels)*	☐	☐	☐	☐	☐	☐	$192.00
Mini Button 1 - Blue *(Velcro)*	☐	☐	☐	☐	☐	☐	$123.00
Mini Button 1 - Green *(Velcro)*	☐	☐	☐	☐	☐	☐	$123.00
Mini Button 1 - Red *(Velcro)*	☐	☐	☐	☐	☐	☐	$123.00
1" Single Switch - Black	☐	☐	☐	☐	☐	☐	$321.00
1" Single Switch - Blue	☐	☐	☐	☐	☐	☐	$321.00
1" Single Switch - Green	☐	☐	☐	☐	☐	☐	$321.00
1" Single Switch - Orange	☐	☐	☐	☐	☐	☐	$321.00
1" Single Switch - Red	☐	☐	☐	☐	☐	☐	$321.00
1" Single Switch - Yellow	☐	☐	☐	☐	☐	☐	$321.00
Electronic Switches							
Flex Switch (Headrest mount only)	☐	☐	☐	☐	☐	☐	$250.000
ProSpot - Proportional Switches							
ProSpot Proportional Switch - Red	☐	☐	☐	☐	☐	☐	$554.00
ProSpot Proportional Switch - Green	☐	☐	☐	☐	☐	☐	$554.00
ProSpot Proportional Switch - Blue	☐	☐	☐	☐	☐	☐	$554.00
ProSpot Proportional Switch - Yellow	☐	☐	☐	☐	☐	☐	$554.00
ProSpot Proportional Switch - Orange	☐	☐	☐	☐	☐	☐	$554.00
Section B – Switch Mounting Position							
Link-It to the Backrest *(Up to 1 per side)*	☐	☐	☐	☐	☐	☐	$759.00
Link-It to the Headrest - Gooseneck *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$759.00
Link-It to the Fixed Joystick Mount *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$578.00
Link-It to the Parallel Swing-Away Joystick Mount *(Up to 2 per side)*	☐	☐	☐	☐	☐	☐	$578.00
Link-It to the Seat Rail *(Up to 2 per side)* *(Not avail w/knee & thigh supports)*	☐	☐	☐	☐	☐	☐	$460.00
Omit Mount *(includes Velcro)*	☐	☐	☐	☐	☐	☐	NC
Section C – Switch Mounting Location							
Right side of chair	☐	☐	☐	☐	☐	☐	NC
Left side of chair	☐	☐	☐	☐	☐	☐	NC

Joystick Integrated Switches *(Pair)*

☐	300JO03	–	Proximity Switches on Flex Mount, Pair *(Used with VersaGuide/VersaGuide EZ)*	$1,161.00

JAY® SYNCRA™ SEATING MEASURING GUIDE

JAY® SYNCRA™ BACK



Backrest Width	16"	18"	20"	22"
Max User Trunk Width	13.5"	15.5"	17.5"	19.5"


JAY® Swing-Away Lateral Thoracic Adjustment Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan
16"	10" to 24"	12" to 19.5"
18"	12" to 26"	12.5" to 21.5"
20"	14" to 28"	12.5" to 22"
22"	16" to 30"	13.5" to 23"

Syncra™ Fixed Lateral Support Accessory Range

Syncra Back Width	Width Range	Height range from bottom of shell to seat pan	Contour Depth Range
16"	16.4" to 16.8"	6" to 8"	3" to 5"
18"	18.4" to 18.8"	6" to 8"	3" to 5"
20"	20.4" to 20.8"	6" to 8"	3" to 5"
22"	22.4" to 22.8"	6" to 8"	3" to 5"

Syncra™ Back Shell Height Diagram



JAY® SYNCRA™ CUSHION



JAY® Lateral Thigh Support Adjustment Range

Syncra Cushion Width	Width Range
16"	8.5" to 16.0"
17"	9.5" to 17.0"
18"	10.5" to 18.0"
19"	11.5" to 19.0"
20"	12.5" to 20.0"
22"	14.5" to 22.0"

Comments for special build

built4me

Q500 M with SEDEO® PRO™	$
Base Price	
Total of additional equipment	
Packaging and transport	
Total Price	

Signature

Send E-Mail